



TW Surgical

Hand & Micro Instrument

Micro Fine Mini Serrefine



TWS 58 - 8010/50	Micro Fine Mini Serrefine, 50mm Length, 10mm Jaw Length 1.2mm Jaw Wide ,Straight, Delicate, Stainless Steel
TWS 58 - 8012/55	Micro Fine Mini Serrefine, 55mm Length, 12mm Jaw Length 1.2mm Jaw Wide ,Straight, Delicate, Stainless Steel
TWS 58 - 8016/60	Micro Fine Mini Serrefine, 60mm Length, 16mm Jaw Length 1.2mm Jaw Wide ,Straight, Delicate, Stainless Steel
TWS 58 - 8020/65	Micro Fine Mini Serrefine, 65mm Length, 20mm Jaw Length 1.2mm Jaw Wide ,Straight, Delicate, Stainless Steel



TWS 58 – 8108/43	Micro Fine Mini Serrefine, 43mm Length 8mm Jaw Light Curved, Extra Narrow and Delicate 1mm jaw wide
TWS 58 – 8112/45	Micro Fine Mini Serrefine, 45mm Length 12mm Jaw Light Curved, Extra Narrow and Delicate 1.2mm jaw wide
TWS 58 – 8116/50	Micro Fine Mini Serrefine, 50mm Length 16mm Jaw Light Curved, Extra Narrow and Delicate 1.2mm jaw wide
TWS 58 – 8120/55	Micro Fine Mini Serrefine, 55mm Length 20mm Jaw Light Curved, Extra Narrow and Delicate 1.2mm jaw wide


Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom

info@twrightmedical.com

Hand & Micro Instrument

Micro Fine Mini Serrefine



TWS 58 - 8063/01	Micro Fine Mini Serrefine, 55mm Length Light Curve Tip, Micro Fine Serrations, Stainless Steel
------------------	---



TWS 58 - 8063/02	Micro Fine Mini Serrefine, 50mm Length Strong Curve Tip, Micro Fine Serrations, Stainless Steel
------------------	--


Hand & Micro Instrument

Micro Fine Mini Serrefine



TWS 58 - 8063/03

Micro Fine Mini Serrefine, 50mm Length
Right Angle Tip 5mm, Micro Fine Serrations, Stainless Steel



TWS 58 - 8063/04

Micro Fine Mini Serrefine, 50mm Length
Right Angle Tip 10mm, Micro Fine Serrations, Stainless Steel