

Hand & Micro Instrument

Micro Needle Holders are used to handle fine needles during suturing procedures.



Plain Jaws, Stainless Steel 0.5mm Wide, 140mm Long

TWS 58 – 8022/05	Micro Needle Holders 140 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Straight Jaw, Straight, Plain, Stainless Steel
TWS 58 – 8023/05	Micro Needle Holders 140 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Curved Jaw, Curved, Plain, Stainless Steel


Diamond Jaws, Stainless Steel 0.5mm Wide, 140mm Long

TWS 58 – 8024/05	Micro Needle Holders 140 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Straight Diamond Jaw, Stainless Steel
TWS 58 – 8025/05	Micro Needle Holders 140 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Curved Diamond Jaw, Stainless Steel

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom
 info@twrightmedical.com

Hand & Micro Instrument




TWS 58 – 8032/05	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8033/05	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Curved Jaw, Straight, Stainless Steel
TWS 58 – 8034/05	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.5mm wide, Angled Jaw, Straight, Stainless Steel

TWS 58 – 8132/075	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.75mm wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8133/075	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.75mm wide, Curved Jaw, Straight, Stainless Steel
TWS 58 – 8134/075	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.75mm wide, Angled Jaw, Straight, Stainless Steel

TWS 58 – 8232/10	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 1.0mm wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8233/10	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 1.0mm wide, Curved Jaw, Straight, Stainless Steel
TWS 58 – 8234/10	Micro Needle Holders 135 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 1.0mm wide, Angled Jaw, Straight, Stainless Steel

Hand & Micro Instrument



Standard, Stainless Steel 0.4mm and 0.8mm Wide, 180mm Long

TWS 58 – 8332/04	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8333/04	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Curved Jaw, Straight, Stainless Steel


TWS 58 – 8432/08	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8433/08	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Curved Jaw, Straight, Stainless Steel

Diamond, Stainless Steel 0.4mm and 0.8mm Wide, 180mm Long

TWS 58 – 8532/04	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Straight Diamond Jaw, Straight, Stainless Steel
TWS 58 – 8533/04	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Curved Diamond Jaw, Straight, Stainless Steel

TWS 58 – 8632/08	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Straight Diamond Jaw, Straight, Stainless Steel
TWS 58 – 8633/08	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Curved Diamond Jaw, Straight, Stainless Steel

Hand & Micro Instrument



Standard, Titanium 0.4mm and 0.8mm Wide, 180mm Long

TWS 58 – 8332/04T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Straight Jaw, Straight, Titanium
TWS 58 – 8333/04T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Curved Jaw, Straight, Titanium


TWS 58 – 8432/08T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Straight Jaw, Straight, Titanium
TWS 58 – 8433/08T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Curved Jaw, Straight, Titanium

Diamond, Titanium 0.4mm and 0.8mm Wide, 180mm Long

TWS 58 – 8532/04T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Straight Diamond Jaw, Straight, Titanium
TWS 58 – 8533/04T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Extra Fine Tapered Tip 0.4mm Wide, Curved Diamond Jaw, Straight, Titanium

TWS 58 – 8632/08T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Straight Diamond Jaw, Straight, Titanium
TWS 58 – 8633/08T	Micro Needle Holders 180 mm, Round Knurled Handle, Tapered Jaws, Fine Tapered Tip 0.8mm Wide, Curved Diamond Jaw, Straight, Titanium

Hand & Micro Instrument





Standard Smooth Carbide Insert, Long Tapered Jaw, Stainless Steel 1.0mm Wide, 180mm Long

TWS 58 – 8732/04	Micro Needle Holders 180 mm, Round Knurled Handle, Long Tapered Jaws, Fine Tapered Tip 1.0mm Wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8733/04	Micro Needle Holders 180 mm, Round Knurled Handle, Long Tapered Jaws, Fine Tapered Tip 1.0mm Wide, Curved Jaw, Straight, Stainless Steel

Standard Serrated Carbide Insert, Long Tapered Jaw, Stainless Steel 1.0mm Wide, 180mm Long

TWS 58 – 8832/08	Micro Needle Holders 180 mm, Round Knurled Handle, Long Tapered Jaws, Fine Tapered Tip 1.0mm Wide, Straight Jaw, Straight, Stainless Steel
TWS 58 – 8833/08	Micro Needle Holders 180 mm, Round Knurled Handle, Long Tapered Jaws, Fine Tapered Tip 1.0mm Wide, Curved Jaw, Straight, Stainless Steel

Micro Needle Holders with Scissors



Plain Jaw, Stainless Steel 0.5mm Wide, 5mm Cutting Length, 140mm Long

TWS 58 - 8030/01	Micro Needle Holders with Scissors 140 mm, Round Knurled Handle, Fine Tapered Tip, 5mm Cutting Length Straight Plain Jaw, Straight, Stainless Steel
TWS 58 – 8030/02	Micro Needle Holders with Scissors 140 mm, Round Knurled Handle, Fine Tapered Tip, 5mm Cutting Length Curved Plain Jaw, Straight, Stainless Steel

Diamond Jaw, Stainless Steel 0.5mm Wide, 5mm Cutting Length, 140mm Long

TWS 58 - 8031/01	Micro Needle Holders with Scissors 140 mm, Round Knurled Handle, Fine Tapered Tip, 5mm Cutting Length Straight Diamond Jaw, Straight, Stainless Steel
TWS 58 – 8031/02	Micro Needle Holders with Scissors 140 mm, Round Knurled Handle, Fine Tapered Tip, 5mm Cutting Length Curved Diamond Jaw, Straight, Stainless Steel