



TW
Surgical

PLASTIC STERILIZATION TRAYS



**UK
CA** **CE**



TW Surgical

Treasure Wright (U.K.) Ltd.



www.twrightmedical.com



info@twrightmedical.com

Sterilization Trays and Products




TW SURGICAL

Plastic Sterilization Trays

- Plastic trays are sterilizable by all standard methods - Autoclave/dry heat up to 360°F, ETO, and Cold Solutions.
- Stackable - Bases feature integral feet and lids feature corner ridges to allow for easier handling in autoclaves and storage.
- All bases have pebble-like surface to prevent mats from sticking.
- Single- and double-level trays are available to handle different instrument quantities. Double Level trays come with a removable upper-level insert and mat.
- Tray dimensions designed to accommodate a variety of instrument sizes and shapes.

Mini Trays

Product Photo	Product Number	Description
	TWS 1300A	1.5" x 2.7" x 0.75" Includes Base, Mat & Lid
	TWS 1310A	1.5" x 2.7" x 1" Includes Deep Base, Mat & Lid
	TWS 1320A	1.5" x 2.7" x 1.25" Includes Deepest Base, Mat & Lid
	TWS 2600A	2.5" x 6" x 0.75" Includes Base, Lid and Mat
	TWS 2610A	2.5" x 6" x 1.25" Includes Base, Lid and Mat



Plastic Sterilization Trays

10-Inch Trays - Featured Products


Enhanced features of the TWS-series trays include a slotted lid and raised node base surface to facilitate the circulation of steam and aid in drying.



TWS 6100A

6" x 10" x 0.75"

Includes Base, Lid and Mat



TWS 6120A

6" x 10" x 1.5" - Double Tray

Includes Base, Insert Tray, and Lid

Sterilization Trays and Products









TW SURGICAL

Plastic Sterilization Trays

Enhanced features of the TWS-series trays include a slotted lid and raised node base surface to facilitate the circulation of steam and aid in drying.

10-Inch Trays

Product Photo	Product Number	Description
	TWS 2700A	2.5" x 7.5" x 0.75" Includes Base, Lid and Mat
	TWS 4600A	4" x 6.5" x 0.75" Includes Base, Lid and Mat
	TWS 4700A	4" x 7.5" x 0.75" Includes Base, Lid and Mat
	TWS 4710A	4" x 7.5" x 1.5" Includes Base, Lid and Mat
	TWS 6124A	6" x 10" x 1.5" Includes Base, Mat and Lid
	TWS 1500A	10" x 15" x 0.75" Includes Base, Lid and Mat



Sterilization Trays and Products









TW SURGICAL

Plastic Sterilization Trays

10-Inch Trays

Product Photo	Product Number	Description
	TWS 1520A	10" x 15" x 1.5" - Complete Set Includes Deep Base, Insert Tray, 2 Mats, and Lid
	TWS 1522A	10" x 15" x 1.5" Includes Deep Base, Base Mat, and Lid

Scope Trays

Product Photo	Product Number	Description
	TWS 6500A	2.6" x 12" x 1.5" Includes Base, Lid (Yellow), and 2 Silicone Bars
	TWS 6510A	2.6" x 16" x 1.5" Includes Base, Lid (Yellow), and 2 Silicone Bars
	TWS 6512A	2.6" x 16" x 1.5" Includes Base, Lid (Yellow), and Mat
	TWS 6520A	3" x 22" x 1.5" Includes Base, Mat & 3 Bars
	TWS 6550A	6.5" x 18" x 1.5" Includes Base, Lid and 6 Bars
	TWS 6552A	6.5" x 18" x 1.5" Includes Base, Lid & Mat




Sterilization Trays and Products




TW SURGICAL

Plastic Sterilization Trays

Scope Trays

Product Photo	Product Number	Description
	TWS 6560A	6.5" x 18" x 3" Includes Base, Lid & Bars
	TWS 6502A	2.6" x 12" x 1.5" Includes Base, Lid, and Mat
	TWS 6562A	6.5" x 18" x 3" Includes Deep Base, Lid and Mat

Silicone Mats

Product Photo	Product Number	Description
	TWS 1380M	13" x 18" Perforated Silicone Mat. Light Blue.

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com








Sterilization Trays and Products



TW SURGICAL

Micro Sterilization Trays

- Lid retainers secure instruments in place and rounded corners increase sterile field integrity.
- The gold-anodized finish minimizes eye fatigue during instrument tip identification.

Product Photo	Product Number	Description
	TWS 2000100007	Exterior: 10.5" x 10.5" x 2.5" Interior: 9.84" x 10.2" x 2.32" Two levels with pin mat and retainers
	TWS 2000100008	Exterior: 10.5" x 10.5" x 1.5" Interior: 9.84" x 10.2" x 1.23" One level with pin mat and retainers
	TWS 2000100086	Exterior: 10.5" x 10.5" x 1.5" Interior: 9.84" x 10.2" x 1.23" Full pin mat version.
	TWS 2000100009	Exterior: 10.5" x 15" x 1.5" Interior: 9.84" x 14.68" x 1.23" One level with pin mat and retainers
	TWS 2000100106	Exterior: 10.5" x 15" x 1.5" Interior: 9.84" x 14.68" x 1.23" Full pin mat version.
	TWS 2000100014	Exterior: 10.5" x 15" x 2.5" Interior: 9.84" x 14.6" x 2.32" Two levels with pin mat, retainers, short tray, and full height accessory area.
	TWS 2000100024	Exterior: 10.5" x 19.5" x 2.5" Interior: 9.84" x 19" x 2.32" Two levels with pin mat, retainers, long tray, and half height accessory area.



Micro Sterilization Trays

Micro Trays - Featured Products



TWS 2000100012

Exterior: 10.5" x 15" x 2.5"
Interior: 9.84" x 14.60" x 2.32"

Two levels with pin mat, retainers, short tray, and full height accessory area.



TWS 2000100013

Exterior: 10.5" x 19.5" x 2.5"
Interior: 9.84" x 19" x 2.32"

Two levels with pin mat, retainers, long tray, and half height accessory area.

Sterilization Trays and Products



TW SURGICAL

Micro Sterilization Trays


The TWS-AST series trays provide a faster drying time than plastic trays. These metal sterilization packaging systems are approximately 400 times more effective than plastic in transferring heat during the sterilization process.

Product Photo	Product Number	Description
	TWS AST2700A	Aluminum Sterilization Tray with Mat - Small 2" x 7" x 0.8"
	TWS AST4700A	Aluminum Sterilization Tray with Mat - Medium 4" x 7.5" x 0.75"
	TWS AST1500A	Aluminum Sterilization Tray with Mat - X-Large 15" x 10" x 0.75"
	TWS AST1520A	Aluminum Sterilization Tray with Mat - X-Large - Double 10" x 15" x 1.5"
	TWS AST61325	Aluminum Sterilization Tray - Double Layer - Deep 6" x 10" x 1.5"
	TWS AST14815A	Aluminum Sterilization Tray with Mats Double layer with side compartment. 6" x 10" x 3.25" (Insert tray is 0.5" high)



Micro Sterilization Trays

Micro Trays - Featured Products



TWS AST6100A

Large Aluminum Sterilization Tray with Mat

6" x 10" x 0.75"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



Micro Sterilization Trays

Micro Trays - Featured Products



TWS AST6120A

Aluminum Sterilization Tray with Mat - Double Layer

6" x 10" x 1.5"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

Sterilization Trays and Products



TW SURGICAL

Scope and Lap Trays

- Nylon coated handle compartment for easy access.
- Moveable silicone inserts can be used to accommodate any size lap instrument (pediatric and adult).
- The top-side handles allow for easy transfer, improving mechanics and reducing potential contamination.

Product Photo	Product Number	Description
	TWS 2000100001	Camera case with brackets Exterior: 9" x 10" x 13" Interior: 8.47" x 8.68" x 2.55"
	TWS 2000100021	Light Cord Case Exterior: 9" x 9" x 1.5" Interior: 8.57" x 8.72" x 1.18"
	TWS 2000100011	Scope Case - 11" - Without Brackets Exterior: 3" x 11" x 2.5" Interior: 2.93" x 10.93" x 2.46"
	TWS 2000100019	Scope Case - 18" - Without Brackets Exterior: 3" x 18" x 2.5" Interior: 2.93" x 17.93" x 2.46"
	TWS 2000100004	Scope Case - Two Levels without brackets / dividers Exterior: 9" x 20" x 4" Interior: 8.47" x 18.68" x 3.55"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com



Scope and Lap Trays

Product Photo	Product Number	Description
	TWS 2000100022	Single Scope Tray with Brackets Exterior: 3.5" x 26" x 3.5" Interior: 3" x 25" x 3"
	TWS 2000100005	Scope Case - Two Levels with brackets / camera and light cord compartments Exterior: 9" x 20" x 4" Interior: 8.47" x 18.68" x 3.55"
	TWS 2000100006	Scope Case with Brackets - Holds 2 to 4 Scopes Exterior: 4" x 17.5" x 1.5" Interior: 3.57" x 17.19" x 1.35"
	TWS 2000100023	Case for 5 Scopes - with Brackets External: 6.5" x 31" x 3" Internal: 6" x 30" x 2.5"
	TWS 2000100035	Surgical Case without Brackets / Pin Mat Exterior: 10" x 21.5" x 2.5" Interior: 9.72" x 20.65" x 2.42"
	TWS 2000100015	Lap Chole Case - Two levels with insert and brackets for instruments up to 35 cm (13.5"). Shown with optional pin mat in the base for additional instrument storage. Exterior: 11" x 23" x 8" Interior: 10.7" x 21.62" x 7.66"






Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

Sterilization Trays and Products



TW SURGICAL

Scope and Lap Trays

Product Photo	Product Number	Description
	TWS 2000100017	Lap Chole Case - Two levels with insert and brackets for instruments up to 65 cm (26"). Shown with optional pin mat in the base for additional instrument storage. Exterior: 11" x 27" x 8" Interior: 10.7" x 25.62" x 7.66"
	TWS 2000100018	Insert for Lap Chole Case - This special assembly holds 14 instruments. It fits inside a 10" x 20" x 6" container. Exterior: 9.5" x 20" x 5.5"
	TWS 2000100026	Insert for Lap Chole Case - Holds 14 instruments. 8" A-Type Container. Fits inside a 10" x 20" x 6" container. Exterior: 9.5" x 20" x 5.5"
	TWS 2000100079	Insert for Lap Chole Case - Holds 14 instruments. Fits inside a 10" x 20" x 6" container. Exterior: 9.9" x 19.5" x 4.9"
	TWS 5200002	Scope Protector 18 centimeters long.

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com




Sterilization Trays and Products



TW SURGICAL

General Instrument Trays

- General surgical instrument trays are available in stainless steel and anodized aluminum and feature medical grade silicone.
- Features include pin mats, insert trays, support brackets, center lid retainers, sturdy latches and handles.

Product Photo	Product Number	Description
	TWS 2000100002	Surgical Utility Case with full pin mat and center lid retainer. Exterior: 5" x 10.5" x 2" Interior: 4.5" x 10.29" x 1.53"
	TWS 2000100003	Surgical Utility Case Exterior: 5" x 10.5" x 2" Interior: 4.5" x 10.29" x 1.53"
	TWS 2000100020	Surgical Case without option tray / support brackets / pin mat Exterior: 11" x 23" x 5" Interior: 10.7" x 21.62" x 3.41"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

GENERAL INSTRUMENT TRAYS

Sterilization Trays and Products



TW SURGICAL

General Instrument Trays

Product Photo	Product Number	Description
	TWS 2000100030	Surgical case without optional tray / support brackets / pin mat Exterior: 11" x 23" x 5" Interior: 10.7" x 21.62" x 4.91"
	TWS 2000100031	Surgical case with optional tray / support brackets but without pin mat Exterior: 11" x 23" x 5" Interior: 10.7" x 21.62" x 4.91"
	TWS 2000100032	Surgical case with two (2) optional trays / support brackets but without pin mat Exterior: 11" x 23" x 5" Interior: 10.7" x 21.62" x 4.91"
	TWS 2000100029	Surgical case with optional tray but without support brackets / pin mat Exterior: 11" x 23" x 3.5" Interior: 10.7" x 21.62" x 3.41"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

GENERAL INSTRUMENT TRAYS

Sterilization Trays and Products



TW SURGICAL

Vacuum and Gravity Containers

Lid and bottom are made from anodized aluminum alloy. Easy to maintain, corrosion resistant, lightweight yet durable, and conducts heat for material drying period.

Gaskets within the retention plates and contain perimter ensure a secure seal of the container, lid and filter.

Snap on filter plates with simple locking mechanisms for quick, secure and easy change of disposable filters.

Round disposable filters with centering hole to assist in proper alignment.

Four (4) ID Tag windows (2 front, 2 back).


Corner design allows for same size container stacking in storage.

ID Tags are available in six colors and can be engraved with up to 18 characters per line.

Low profile latch easily and securely closes the lid over the container bottom.

Two (2) Indicator Card windows (1 front, 1 back). Indicator card has visible bar for clear indication that steam or gas sterilization has been achieved.

Tamper evident lock - identifies sterility breach.



Container Sizes

Colored lids, wire and perforated baskets are available separately for each size vacuum container except for wide body containers (silver lids and wire baskets only).

Half Size 11.25" x 11.75"	Available in these heights: 4 1/8" 5 1/4" 6 7/8" 7 3/4" 10 1/4"
Three Quarter Size 18.25" x 11.25"	Available in these heights: 4 1/8" 5 1/4" 6"
Full Size 23.25" x 11.25"	Available in these heights: 4 1/8" 5 1/4" 5 7/8" 7 3/4" 10 1/4"
Wide Body 23.25" x 15.75"	Available in these heights: 4 3/4" 7"



Wire and Perforated Baskets



Wire Basket and Perforated Basket Sizes

Wire baskets come with “screw on” feet. Perforated baskets come with “snap on” feet. Partition sheets, holding pins and other accessories are also available.

Half Size 11.25" x 11.75"	Available in these heights: 2" 3" 4"
Three Quarter Size 18.25" x 11.25"	Available in these heights: 2" 3" 4"
Full Size 23.25" x 11.25"	Available in these heights: 2" 3" 4"
Wide Body 23.25" x 15.75"	Available in these heights: 2" 4"

Treasure Wright (U.K.) Ltd.
128-160 City Road, London, EC1V 2NX, United Kingdom
info@twrightmedical.com

Our experienced team of surgical instrument specialists help thousands of facilities each year save time and money on instruments, sets, sterilization trays, and more.

Treasure Wright (U.K.) Ltd.

128- 160 City Road, London, EC1V 2NX, United Kingdom.

INFO@TWRIGHTMEDICAL.COM

www.twrightmedical.com