




TW Surgical



# *Sterile Container Systems Product Catalog*





TW Surgical



## TABLE OF CONTENT

---

04 - 07 INTRODUCTION PAGES  
13 - 25 STANDARD MODEL  
26 - 33 PREMIUM MODEL  
34 - 39 VALVE BARRIER MODEL  
40 - 53 AIO MODEL  
54 - 58 PLASMA MODEL  
60 - 64 REMOVAL MODEL  
66 - 71 DENTAL MODEL  
72 - 77 MINI MODEL  
78 - 85 BASKETS  
86 - 97 ACCESSORIES





TW Surgical



## ABOUT US

TW Surgical is a trusted provider of orthopedic implants, sterilization equipment, surgical instruments, and hospital supplies, all crafted to meet international standards. We are committed to innovation, safety, and reliability, supporting healthcare professionals in delivering exceptional patient care. With a global outlook, our mission is to empower medical institutions with solutions that make a real difference.




TW Surgical


## OUR MISSION

TW Surgical specializes in orthopedic, sterilization, and hospital products, offering high-quality solutions designed to meet modern healthcare needs. Our product range includes orthopedic implants, sterilization equipment, surgical instruments, and essential hospital supplies, all manufactured with a focus on safety, innovation, and international standards. With a commitment to reliability and customer satisfaction, TW Surgical supports medical professionals and institutions in delivering better patient care every day.






*BOTTOM  
NON - PERFORATED*




*BOTTOM  
PERFORATED*






TW Surgical





TW Surgical





**Cost-Efficiency:** Although disposable products might seem economical, the need for frequent disposal makes them less cost-effective in the long run. Using sterilization containers proves to be more economical over time. Additionally, disposable wraps used in pouching can be damaged during sterilization or transport, compromising sterility and often requiring a repeat sterilization cycle. Aluminum sterilization containers help prevent these additional costs and save time by providing a more reliable and durable solution.



**Eco-Friendly:** Disposable products generate significant waste due to their paper and plastic content. By contrast, sterilization container systems reduce waste, helping to prevent the environmental pollution associated with disposable materials.



**Portability:** Equipped with stainless steel handles and optional high-temperature-resistant, antibacterial silicone covers, containers can be transported more easily and safely.



**Identification:** Barcodes, QR codes, or Data Matrix codes can be engraved on colored aluminum labels on both sides of the container for easy identification of instrument sets. Additionally, different lid color options allow for clear separation and identification of sets.




**Safety:** Sterilization with container systems offers greater safety than other methods. When stored properly and unopened, sterilization containers can maintain instrument sterility for up to six months.

## THE ADVANTAGES OF USING STERILE CONTAINER SYSTEMS

Sterilization containers offer an efficient alternative to individually packaging surgical instruments with sterilization rolls or disposable wraps. This container-based sterilization method provides several advantages over these traditional approaches.



# ASSEMBLY COMPONENTS





TW Surgical

# 9 COLOR OPTIONS





TW Surgical

# STANDART MODEL STERILIZATION CONTAINERS



# STANDART SERIES

## 1/2 HALF SIZE

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 14-15 are presented as a set including the base and cover.



TW Surgical



### BOTTOM NON-PERFORATED



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
<b>Gray</b>	PC.A 10.10	PC.A 10.13	PC.A 10.15	PC.A 10.20	PC.A 10.26
<b>Yellow</b>	PC.A 20.10	PC.A 20.13	PC.A 20.15	PC.A 20.20	PC.A 20.26
<b>Red</b>	PC.A 30.10	PC.A 30.13	PC.A 30.15	PC.A 30.20	PC.A 30.26
<b>Blue</b>	PC.A 40.10	PC.A 40.13	PC.A 40.15	PC.A 40.20	PC.A 40.26
<b>Green</b>	PC.A 50.10	PC.A 50.13	PC.A 50.15	PC.A 50.20	PC.A 50.26
<b>Black</b>	PC.A 60.10	PC.A 60.13	PC.A 60.15	PC.A 60.20	PC.A 60.26
<b>Orange</b>	PC.A 70.10	PC.A 70.13	PC.A 70.15	PC.A 70.20	PC.A 70.26
<b>Purple</b>	PC.A 80.10	PC.A 80.13	PC.A 80.15	PC.A 80.20	PC.A 80.26
<b>Bronze</b>	PC.A 90.10	PC.A 90.13	PC.A 90.15	PC.A 90.20	PC.A 90.26

# STANDART SERIES

## 1/2 HALF SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 14-15 are presented as a set including the base and cover.



### BOTTOM PERFORATED



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.B 10.10	PC.B 10.13	PC.B 10.15	PC.B 10.20	PC.B 10.26
Yellow	PC.B 20.10	PC.B 20.13	PC.B 20.15	PC.B 20.20	PC.B 20.26
Red	PC.B 30.10	PC.B 30.13	PC.B 30.15	PC.B 30.20	PC.B 30.26
Blue	PC.B 40.10	PC.B 40.13	PC.B 40.15	PC.B 40.20	PC.B 40.26
Green	PC.B 50.10	PC.B 50.13	PC.B 50.15	PC.B 50.20	PC.B 50.26
Black	PC.B 60.10	PC.B 60.13	PC.B 60.15	PC.B 60.20	PC.B 60.26
Orange	PC.B 70.10	PC.B 70.13	PC.B 70.15	PC.B 70.20	PC.B 70.26
Purple	PC.B 80.10	PC.B 80.13	PC.B 80.15	PC.B 80.20	PC.B 80.26
Bronze	PC.B 90.10	PC.B 90.13	PC.B 90.15	PC.B 90.20	PC.B 90.26

# STANDART SERIES

## 1/2 HALF SIZE




- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!


- The article numbers for products on pages 15-16 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



### BOTTOM PERFORATED



# STANDART SERIES

## 1/2 HALF SIZE

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 16-17 are presented as a set including the base and cover.



### STANDART LIDS FOR HALF SIZES

Colours	285 x 280mm
Gray	PC.A 10.01
Yellow	PC.A 20.01
Red	PC.A 30.01
Blue	PC.A 40.01
Green	PC.A 50.01
Black	PC.A 60.01
Orange	PC.A 70.01
Purple	PC.A 80.01
Bronze	PC.A 90.01



### SAFETY LIDS FOR HALF SIZES

Colours	310 x 285mm
Gray	PC.A 10.02
Yellow	PC.A 20.02
Red	PC.A 30.02
Blue	PC.A 40.02
Green	PC.A 50.02
Black	PC.A 60.02
Orange	PC.A 70.02
Purple	PC.A 80.02
Bronze	PC.A 90.02

# STANDART SERIES

## 3/4 MIDDLE SIZE

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 18-19 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.C 10.10	PC.C 10.13	PC.C 10.15
Yellow	PC.C 20.10	PC.C 20.13	PC.C 20.15
Red	PC.C 30.10	PC.C 30.13	PC.C 30.15
Blue	PC.C 40.10	PC.C 40.13	PC.C 40.15
Green	PC.C 50.10	PC.C 50.13	PC.C 50.15
Black	PC.C 60.10	PC.C 60.13	PC.C 60.15
Orange	PC.C 70.10	PC.C 70.13	PC.C 70.15
Purple	PC.C 80.10	PC.C 80.13	PC.C 80.15
Bronze	PC.C 90.10	PC.C 90.13	PC.C 90.15

# STANDART SERIES

## 3/4 MIDDLE SIZE

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 18-19 are presented as a set including the base and cover.



### BOTTOM PERFORATED



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.D 10.10	PC.D 10.13	PC.D 10.15
Yellow	PC.D 20.10	PC.D 20.13	PC.D 20.15
Red	PC.D 30.10	PC.D 30.13	PC.D 30.15
Blue	PC.D 40.10	PC.D 40.13	PC.D 40.15
Green	PC.D 50.10	PC.D 50.13	PC.D 50.15
Black	PC.D 60.10	PC.D 60.13	PC.D 60.15
Orange	PC.D 70.10	PC.D 70.13	PC.D 70.15
Purple	PC.D 80.10	PC.D 80.13	PC.D 80.15
Bronze	PC.D 90.10	PC.D 90.13	PC.D 90.15

# STANDART SERIES

## 3/4 MIDDLE SIZE




- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 20-21 are presented seperated bottom and lid.



### BOTTOM NON-PERFORATED



465 x 280 x 96mm

PC.C 00.10



465 x 280 x 119mm

PC.C 00.13



465 x 280 x 138mm

PC.C 00.15



### BOTTOM PERFORATED



285 x 280 x 96mm

PC.D 00.10



285 x 280 x 119mm

PC.D 00.13



285 x 280 x 138mm

PC.D 00.15



# STANDART SERIES

## 3/4 MIDDLE SIZE

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 20-21 are presented seperated bottom and lid.



### STANDART LIDS FOR MIDDLE SIZES

Colours 465 x 280mm

Gray PC.C 10.01

Yellow PC.C 20.01

Red PC.C 30.01

Blue PC.C 40.01

Green PC.C 50.01

Black PC.C 60.01

Orange PC.C 70.01

Purple PC.C 80.01

Bronze PC.C 90.01



### SAFETY LIDS FOR MIDDLE SIZES

Colours 485 x 285mm

Gray PC.C 10.02

Yellow PC.C 20.02

Red PC.C 30.02

Blue PC.C 40.02

Green PC.C 50.02

Black PC.C 60.02

Orange PC.C 70.02

Purple PC.C 80.02

Bronze PC.C 90.02

# STANDART SERIES

## 1/1 FULL SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 22-23 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
<b>Gray</b>	PC.E 10.10	PC.E 10.13	PC.E 10.15	PC.E 10.20	PC.E 10.26
<b>Yellow</b>	PC.E 20.10	PC.E 20.13	PC.E 20.15	PC.E 20.20	PC.E 20.26
<b>Red</b>	PC.E 30.10	PC.E 30.13	PC.E 30.15	PC.E 30.20	PC.E 30.26
<b>Blue</b>	PC.E 40.10	PC.E 40.13	PC.E 40.15	PC.E 40.20	PC.E 40.26
<b>Green</b>	PC.E 50.10	PC.E 50.13	PC.E 50.15	PC.E 50.20	PC.E 50.26
<b>Black</b>	PC.E 60.10	PC.E 60.13	PC.E 60.15	PC.E 60.20	PC.E 60.26
<b>Orange</b>	PC.E 70.10	PC.E 70.13	PC.E 70.15	PC.E 70.20	PC.E 70.26
<b>Purple</b>	PC.E 80.10	PC.E 80.13	PC.E 80.15	PC.E 80.20	PC.E 80.26
<b>Bronze</b>	PC.E 90.10	PC.E 90.13	PC.E 90.15	PC.E 90.20	PC.E 90.26

# STANDART SERIES

## 1/1 FULL SIZE



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 22-23 are presented as a set including the base and cover.



### BOTTOM PERFORATED



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.F 10.10	PC.F 10.13	PC.F 10.15	PC.F 10.20	PC.F 10.26
Yellow	PC.F 20.10	PC.F 20.13	PC.F 20.15	PC.F 20.20	PC.F 20.26
Red	PC.F 30.10	PC.F 30.13	PC.F 30.15	PC.F 30.20	PC.F 30.26
Blue	PC.F 40.10	PC.F 40.13	PC.F 40.15	PC.F 40.20	PC.F 40.26
Green	PC.F 50.10	PC.F 50.13	PC.F 50.15	PC.F 50.20	PC.F 50.26
Black	PC.F 60.10	PC.F 60.13	PC.F 60.15	PC.F 60.20	PC.F 60.26
Orange	PC.F 70.10	PC.F 70.13	PC.F 70.15	PC.F 70.20	PC.F 70.26
Purple	PC.F 80.10	PC.F 80.13	PC.F 80.15	PC.F 80.20	PC.F 80.26
Bronze	PC.F 90.10	PC.F 90.13	PC.F 90.15	PC.F 90.20	PC.F 90.26

# STANDART SERIES

## 1/1 FULL SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 24-25 are presented seperated bottom and lid.



### BOTTOM NON-PERFORATED



580 x 280 x 96mm

580 x 280 x 119mm

580 x 280 x 138mm

580 x 280 x 185mm

580 x 280 x 245mm

PC.E 00.10

PC.E 00.13

PC.E 00.15

PC.E 00.20

PC.E 00.26

### BOTTOM PERFORATED



580 x 280 x 96mm

580 x 280 x 119mm

580 x 280 x 138mm

580 x 280 x 185mm

580 x 280 x 245mm

PC.F 00.10

PC.F 00.13

PC.F 00.15

PC.F 00.20

PC.F 00.26

# STANDART SERIES

## 1/1 FULL SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.


\*\*\* Attention!

- The article numbers for products on pages 24-25 are presented seperated bottom and lid.



### STANDART LIDS FOR FULL SIZES

Colours	580 x 280mm
Gray	PC.E 10.01
Yellow	PC.E 20.01
Red	PC.E 30.01
Blue	PC.E 40.01
Green	PC.E 50.01
Black	PC.E 60.01
Orange	PC.E 70.01
Purple	PC.E 80.01
Bronze	PC.E 90.01



### SAFETY LIDS FOR FULL SIZES

Colours	605 x 285mm
Gray	PC.E 10.02
Yellow	PC.E 20.02
Red	PC.E 30.02
Blue	PC.E 40.02
Green	PC.E 50.02
Black	PC.E 60.02
Orange	PC.E 70.02
Purple	PC.E 80.02
Bronze	PC.E 90.02




TW Surgical





TW Surgical

# PREMIUM MODEL STERILIZATION CONTAINERS



# PREMIUM SERIES

## 1/2 HALF SIZE



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 28-29 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
<b>Gray</b>	PC.AP 10.10	PC.AP 10.13	PC.AP 10.15	PC.AP 10.20	PC.AP 10.26
<b>Yellow</b>	PC.AP 20.10	PC.AP 20.13	PC.AP 20.15	PC.AP 20.20	PC.AP 20.26
<b>Red</b>	PC.AP 30.10	PC.AP 30.13	PC.AP 30.15	PC.AP 30.20	PC.AP 30.26
<b>Blue</b>	PC.AP 40.10	PC.AP 40.13	PC.AP 40.15	PC.AP 40.20	PC.AP 40.26
<b>Green</b>	PC.AP 50.10	PC.AP 50.13	PC.AP 50.15	PC.AP 50.20	PC.AP 50.26
<b>Black</b>	PC.AP 60.10	PC.AP 60.13	PC.AP 60.15	PC.AP 60.20	PC.AP 60.26
<b>Orange</b>	PC.AP 70.10	PC.AP 70.13	PC.AP 70.15	PC.AP 70.20	PC.AP 70.26
<b>Purple</b>	PC.AP 80.10	PC.AP 80.13	PC.AP 80.15	PC.AP 80.20	PC.AP 80.26
<b>Bronze</b>	PC.AP 90.10	PC.AP 90.13	PC.AP 90.15	PC.AP 90.20	PC.AP 90.26

# PREMIUM SERIES

## 1/2 HALF SIZE



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/2 Half Size.
  - 9 Different colour options available.
  - Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.
- \*\*\* Attention!
- The article numbers for products on pages 28-29 are presented as a set including the base and cover.



### BOTTOM PERFORATED



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.BP 10.10	PC.BP 10.13	PC.BP 10.15	PC.BP 10.20	PC.BP 10.26
Yellow	PC.BP 20.10	PC.BP 20.13	PC.BP 20.15	PC.BP 20.20	PC.BP 20.26
Red	PC.BP 30.10	PC.BP 30.13	PC.BP 30.15	PC.BP 30.20	PC.BP 30.26
Blue	PC.BP 40.10	PC.BP 40.13	PC.BP 40.15	PC.BP 40.20	PC.BP 40.26
Green	PC.BP 50.10	PC.BP 50.13	PC.BP 50.15	PC.BP 50.20	PC.BP 50.26
Black	PC.BP 60.10	PC.BP 60.13	PC.BP 60.15	PC.BP 60.20	PC.BP 60.26
Orange	PC.BP 70.10	PC.BP 70.13	PC.BP 70.15	PC.BP 70.20	PC.BP 70.26
Purple	PC.BP 80.10	PC.BP 80.13	PC.BP 80.15	PC.BP 80.20	PC.BP 80.26
Bronze	PC.BP 90.10	PC.BP 90.13	PC.BP 90.15	PC.BP 90.20	PC.BP 90.26

# PREMIUM SERIES

## 3/4 MIDDLE SIZE



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 30-31 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.CP 10.10	PC.CP 10.13	PC.CP 10.15
Yellow	PC.CP 20.10	PC.CP 20.13	PC.CP 20.15
Red	PC.CP 30.10	PC.CP 30.13	PC.CP 30.15
Blue	PC.CP 40.10	PC.CP 40.13	PC.CP 40.15
Green	PC.CP 50.10	PC.CP 50.13	PC.CP 50.15
Black	PC.CP 60.10	PC.CP 60.13	PC.CP 60.15
Orange	PC.CP 70.10	PC.CP 70.13	PC.CP 70.15
Purple	PC.CP 80.10	PC.CP 80.13	PC.CP 80.15
Bronze	PC.CP 90.10	PC.CP 90.13	PC.CP 90.15

# PREMIUM SERIES

## 3/4 MIDDLE SIZE



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 3/4 Middle Size.
  - 9 Different colour options available.
  - Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.
- \*\*\* Attention!
- The article numbers for products on pages 30-31 are presented as a set including the base and cover.



### BOTTOM PERFORATED



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.DP 10.10	PC.DP 10.13	PC.DP 10.15
Yellow	PC.DP 20.10	PC.DP 20.13	PC.DP 20.15
Red	PC.DP 30.10	PC.DP 30.13	PC.DP 30.15
Blue	PC.DP 40.10	PC.DP 40.13	PC.DP 40.15
Green	PC.DP 50.10	PC.DP 50.13	PC.DP 50.15
Black	PC.DP 60.10	PC.DP 60.13	PC.DP 60.15
Orange	PC.DP 70.10	PC.DP 70.13	PC.DP 70.15
Purple	PC.DP 80.10	PC.DP 80.13	PC.DP 80.15
Bronze	PC.DP 90.10	PC.DP 90.13	PC.DP 90.15

# PREMIUM SERIES

## 1/1 FULL SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 32-33 are presented as a set including the base and cover.



### BOTTOM NON-PERFORATED



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.EP 10.10	PC.EP 10.13	PC.EP 10.15	PC.EP 10.20	PC.EP 10.26
Yellow	PC.EP 20.10	PC.EP 20.13	PC.EP 20.15	PC.EP 20.20	PC.EP 20.26
Red	PC.EP 30.10	PC.EP 30.13	PC.EP 30.15	PC.EP 30.20	PC.EP 30.26
Blue	PC.EP 40.10	PC.EP 40.13	PC.EP 40.15	PC.EP 40.20	PC.EP 40.26
Green	PC.EP 50.10	PC.EP 50.13	PC.EP 50.15	PC.EP 50.20	PC.EP 50.26
Black	PC.EP 60.10	PC.EP 60.13	PC.EP 60.15	PC.EP 60.20	PC.EP 60.26
Orange	PC.EP 70.10	PC.EP 70.13	PC.EP 70.15	PC.EP 70.20	PC.EP 70.26
Purple	PC.EP 80.10	PC.EP 80.13	PC.EP 80.15	PC.EP 80.20	PC.EP 80.26
Bronze	PC.E 90.10	PC.EP 90.13	PC.EP 90.15	PC.EP 90.20	PC.EP 90.26

# PREMIUM SERIES

## 1/1 FULL SIZE




TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of 1/1 Full Size.
- 9 Different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.

\*\*\* Attention!

- The article numbers for products on pages 32-33 are presented as a set including the base and cover.




### BOTTOM PERFORATED



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.FP 10.10	PC.FP 10.13	PC.FP 10.15	PC.FP 10.20	PC.FP 10.26
Yellow	PC.FP 20.10	PC.FP 20.13	PC.FP 20.15	PC.FP 20.20	PC.FP 20.26
Red	PC.FP 30.10	PC.FP 30.13	PC.FP 30.15	PC.FP 30.20	PC.FP 30.26
Blue	PC.FP 40.10	PC.FP 40.13	PC.FP 40.15	PC.FP 40.20	PC.FP 40.26
Green	PC.FP 50.10	PC.FP 50.13	PC.FP 50.15	PC.FP 50.20	PC.FP50.26
Black	PC.FP 60.10	PC.FP 60.13	PC.FP 60.15	PC.FP 60.20	PC.FP 60.26
Orange	PC.FP 70.10	PC.FP 70.13	PC.FP 70.15	PC.FP 70.20	PC.FP 70.26
Purple	PC.FP 80.10	PC.FP 80.13	PC.FP 80.15	PC.FP 80.20	PC.FP 80.26
Bronze	PC.FP 90.10	PC.FP 90.13	PC.FP 90.15	PC.FP 90.20	PC.FP 90.26



TW Surgical






TW Surgical


# CLASSIC VALVE BIOBARRIER MODEL STERILIZATION CONTAINERS




# CLASSIC VALVE BIOBARRIER SERIES STERILIZATION CONTAINERS




Valves are in the closed position before sterilization, as shown in the figure.




Valves are automatically opened by the steam pressure in the steam sterilizer and the sterilization cycle begins.



When the steam pressure in the steam sterilizer returns to normal, the valves are opened in the opposite direction and the container internal pressure and sterilizer pressure are equalized.



Valves become completely closed and air circulation is blocked after the sterilization cycle.




# CLASSIC VALVE BIOBARRIER SERIES

## 1/2 HALF SIZE



- Valve Barrier series containers are working without filter.
- 9 Different colour options available.



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.AB 10.10	PC.AB 10.13	PC.AB 10.15	PC.AB 10.20	PC.AB 10.26
Yellow	PC.AB 20.10	PC.AB 20.13	PC.AB 20.15	PC.AB 20.20	PC.AB 20.26
Red	PC.AB 30.10	PC.AB 30.13	PC.AB 30.15	PC.AB 30.20	PC.AB 30.26
Blue	PC.AB 40.10	PC.AB 40.13	PC.AB 40.15	PC.AB 40.20	PC.AB 40.26
Green	PC.AB 50.10	PC.AB 50.13	PC.AB 50.15	PC.AB 50.20	PC.AB 50.26
Black	PC.AB 60.10	PC.AB 60.13	PC.AB 60.15	PC.AB 60.20	PC.AB 60.26
Orange	PC.AB 70.10	PC.AB 70.13	PC.AB 70.15	PC.AB 70.20	PC.AB 70.26
Purple	PC.AB 80.10	PC.AB 80.13	PC.AB 80.15	PC.AB 80.20	PC.AB 80.26
Bronze	PC.AB 90.10	PC.AB 90.13	PC.AB 90.15	PC.AB 90.20	PC.AB 90.26

# CLASSIC VALVE BIOBARRIER SERIES

## 3/4 MIDDLE SIZE

- Valve Barrier series containers are working without filter.
- 9 Different colour options available.



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.BB 10.10	PC.BB 10.13	PC.BB 10.15
Yellow	PC.BB 20.10	PC.BB 20.13	PC.BB 20.15
Red	PC.BB 30.10	PC.BB 30.13	PC.BB 30.15
Blue	PC.BB 40.10	PC.BB 40.13	PC.BB 40.15
Green	PC.BB 50.10	PC.BB 50.13	PC.BB 50.15
Black	PC.BB 60.10	PC.BB 60.13	PC.BB 60.15
Orange	PC.BB 70.10	PC.BB 70.13	PC.BB 70.15
Purple	PC.BB 80.10	PC.BB 80.13	PC.BB 80.15
Bronze	PC.BB 90.10	PC.BB 90.13	PC.BB 90.15


# CLASSIC VALVE BIOBARRIER SERIES

## 1/1 FULL SIZE



TW Surgical

- Valve Barrier series containers are working without filter.
- 9 Different colour options available.



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.CB 10.10	PC.CB 10.13	PC.CB 10.15	PC.CB 10.20	PC.CB 10.26
Yellow	PC.CB 20.10	PC.CB 20.13	PC.CB 20.15	PC.CB 20.20	PC.CB 20.26
Red	PC.CB 30.10	PC.CB 30.13	PC.CB 30.15	PC.CB 30.20	PC.CB 30.26
Blue	PC.CB 40.10	PC.CB 40.13	PC.CB 40.15	PC.CB 40.20	PC.CB 40.26
Green	PC.CB 50.10	PC.CB 50.13	PC.CB 50.15	PC.CB 50.20	PC.CB 50.26
Black	PC.CB 60.10	PC.CB 60.13	PC.CB 60.15	PC.CB 60.20	PC.CB 60.26
Orange	PC.CB 70.10	PC.CB 70.13	PC.CB 70.15	PC.CB 70.20	PC.CB 70.26
Purple	PC.CB 80.10	PC.CB 80.13	PC.CB 80.15	PC.CB 80.20	PC.CB 80.26
Bronze	PC.CB 90.10	PC.CB 90.13	PC.CB 90.15	PC.CB 90.20	PC.CB 90.26



TW Surgical





TW Surgical

# SMART AIO MODEL STERILIZATION CONTAINERS



# SMART AIO SERIES VALVE OPTION





# SMART AIO SERIES - Valve Option

## 1/2 HALF SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
<b>Gray</b>	PC.AV 10.10	PC.AV 10.13	PC.AV 10.15	PC.AV 10.20	PC.AV 10.26
<b>Yellow</b>	PC.AV 20.10	PC.AV 20.13	PC.AV 20.15	PC.AV 20.20	PC.AV 20.26
<b>Red</b>	PC.AV 30.10	PC.AV 30.13	PC.AV 30.15	PC.AV 30.20	PC.AV 30.26
<b>Blue</b>	PC.AV 40.10	PC.AV 40.13	PC.AV 40.15	PC.AV 40.20	PC.AV 40.26
<b>Green</b>	PC.AV 50.10	PC.AV 50.13	PC.AV 50.15	PC.AV 50.20	PC.AV 50.26
<b>Black</b>	PC.AV 60.10	PC.AV 60.13	PC.AV 60.15	PC.AV 60.20	PC.AV 60.26
<b>Orange</b>	PC.AV 70.10	PC.AV 70.13	PC.AV 70.15	PC.AV 70.20	PC.AV 70.26
<b>Purple</b>	PC.AV 80.10	PC.AV 80.13	PC.AV 80.15	PC.AV 80.20	PC.AV 80.26
<b>Bronze</b>	PC.AV 90.10	PC.AV 90.13	PC.AV 90.15	PC.AV 90.20	PC.AV 90.26

# SMART AIO SERIES - Valve Option

## 3/4 MIDDLE SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.




Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.BV 10.10	PC.BV 10.13	PC.BV 10.15
Yellow	PC.BV 20.10	PC.BV 20.13	PC.BV 20.15
Red	PC.BV 30.10	PC.BV 30.13	PC.BV 30.15
Blue	PC.BV 40.10	PC.BV 40.13	PC.BV 40.15
Green	PC.BV 50.10	PC.BV 50.13	PC.BV 50.15
Black	PC.BV 60.10	PC.BV 60.13	PC.BV 60.15
Orange	PC.BV 70.10	PC.BV 70.13	PC.BV 70.15
Purple	PC.BV 80.10	PC.BV 80.13	PC.BV 80.15
Bronze	PC.BV 90.10	PC.BV 90.13	PC.BV 90.15

# SMART AIO SERIES - Valve Option

## 1/1 FULL SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.CV 10.10	PC.CV 10.13	PC.CV 10.15	PC.CV 10.20	PC.CV 10.26
Yellow	PC.CV 20.10	PC.CV 20.13	PC.CV 20.15	PC.CV 20.20	PC.CV 20.26
Red	PC.CV 30.10	PC.CV 30.13	PC.CV 30.15	PC.CV 30.20	PC.CV 30.26
Blue	PC.CV 40.10	PC.CV 40.13	PC.CV 40.15	PC.CV 40.20	PC.CV 40.26
Green	PC.CV 50.10	PC.CV 50.13	PC.CV 50.15	PC.CV 50.20	PC.CV 50.26
Black	PC.CV 60.10	PC.CV 60.13	PC.CV 60.15	PC.CV 60.20	PC.CV 60.26
Orange	PC.CV 70.10	PC.CV 70.13	PC.CV 70.15	PC.CV 70.20	PC.CV 70.26
Purple	PC.CV 80.10	PC.CV 80.13	PC.CV 80.15	PC.CV 80.20	PC.CV 80.26
Bronze	PC.CV 90.10	PC.CV 90.13	PC.CV 90.15	PC.CV 90.20	PC.CV 90.26

**AIO SERIES**

# **CASSETTE FILTER OPTION**





TW Surgical

# AIO SERIES - Casette Filter Option

## 1/2 HALF SIZE



- AIO Series - Casette Option containers are working with cassette filters.
- AIO series containers filter systems can be detachable, this containers can work with PTFE, Textile and Paper filter types.
- 9 Different colour options available.



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.AC 10.10	PC.AC 10.13	PC.AC 10.15	PC.AC 10.20	PC.AC 10.26
Yellow	PC.AC 20.10	PC.AC 20.13	PC.AC 20.15	PC.AC 20.20	PC.AC 20.26
Red	PC.AC 30.10	PC.AC 30.13	PC.AC 30.15	PC.AC 30.20	PC.AC 30.26
Blue	PC.AC 40.10	PC.AC 40.13	PC.AC 40.15	PC.AC 40.20	PC.AC 40.26
Green	PC.AC 50.10	PC.AC 50.13	PC.AC 50.15	PC.AC 50.20	PC.AC 50.26
Black	PC.AC 60.10	PC.AC 60.13	PC.AC 60.15	PC.AC 60.20	PC.AC 60.26
Orange	PC.AC 70.10	PC.AC 70.13	PC.AC 70.15	PC.AC 70.20	PC.AC 70.26
Purple	PC.AC 80.10	PC.AC 80.13	PC.AC 80.15	PC.AC 80.20	PC.AC 80.26
Bronze	PC.AC 90.10	PC.AC 90.13	PC.AC 90.15	PC.AC 90.20	PC.AC 90.26

# AIO SERIES - Casette Filter Option

## 3/4 MIDDLE SIZE



- AIO Series - Casette Option containers are working with cassette filters.
- AIO series containers filter systems can be detachable, this containers can work with PTFE, Textile and Paper filter types.
- 9 Different colour options available.





Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.BC 10.10	PC.BC 10.13	PC.BC 10.15
Yellow	PC.BC 20.10	PC.BC 20.13	PC.BC 20.15
Red	PC.BC 30.10	PC.BC 30.13	PC.BC 30.15
Blue	PC.BC 40.10	PC.BC 40.13	PC.BC 40.15
Green	PC.BC 50.10	PC.BC 50.13	PC.BC 50.15
Black	PC.BC 60.10	PC.BC 60.13	PC.BC 60.15
Orange	PC.BC 70.10	PC.BC 70.13	PC.BC 70.15
Purple	PC.BC 80.10	PC.BC 80.13	PC.BC 80.15
Bronze	PC.BC 90.10	PC.BC 90.13	PC.BC 90.15

# AIO SERIES - Casette Filter Option

## 1/1 FULL SIZE



- AIO Series - Casette Option containers are working with cassette filters.
- AIO series containers filter systems can be detachable, this containers can work with PTFE, Textile and Paper filter types.
- 9 Different colour options available.




Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.CC 10.10	PC.CC 10.13	PC.CC 10.15	PC.CC 10.20	PC.CC 10.26
Yellow	PC.CC 20.10	PC.CC 20.13	PC.CC 20.15	PC.CC 20.20	PC.CC 20.26
Red	PC.CC 30.10	PC.CC 30.13	PC.CC 30.15	PC.CC 30.20	PC.CC 30.26
Blue	PC.CC 40.10	PC.CC 40.13	PC.CC 40.15	PC.CC 40.20	PC.CC 40.26
Green	PC.CC 50.10	PC.CC 50.13	PC.CC 50.15	PC.CC 50.20	PC.CC 50.26
Black	PC.CC 60.10	PC.CC 60.13	PC.CC 60.15	PC.CC 60.20	PC.CC 60.26
Orange	PC.CC 70.10	PC.CC 70.13	PC.CC 70.15	PC.CC 70.20	PC.CC 70.26
Purple	PC.CC 80.10	PC.CC 80.13	PC.CC 80.15	PC.CC 80.20	PC.CC 80.26
Bronze	PC.CC 90.10	PC.CC 90.13	PC.CC 90.15	PC.CC 90.20	PC.CC 90.26

**AIO SERIES**

# **VALVE BARRIER OPTION**



TW Surgical





# AIO SERIES - Valve Barrier Option

## 1/2 HALF SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.AS 10.10	PC.AS 10.13	PC.AS 10.15	PC.AS 10.20	PC.AS 10.26
Yellow	PC.AS 20.10	PC.AS 20.13	PC.AS 20.15	PC.AS 20.20	PC.AS 20.26
Red	PC.AS 30.10	PC.AS 30.13	PC.AS 30.15	PC.AS 30.20	PC.AS 30.26
Blue	PC.AS 40.10	PC.AS 40.13	PC.AS 40.15	PC.AS 40.20	PC.AS 40.26
Green	PC.AS 50.10	PC.AS 50.13	PC.AS 50.15	PC.AS 50.20	PC.AS 50.26
Black	PC.AS 60.10	PC.AS 60.13	PC.AS 60.15	PC.AS 60.20	PC.AS 60.26
Orange	PC.AS 70.10	PC.AS 70.13	PC.AS 70.15	PC.AS 70.20	PC.AS 70.26
Purple	PC.AS 80.10	PC.AS 80.13	PC.AS 80.15	PC.AS 80.20	PC.AS 80.26
Bronze	PC.AS 90.10	PC.AC 90.13	PC.AS 90.15	PC.AS 90.20	PC.AS 90.26

# AIO SERIES - Valve Barrier Option

## 3/4 MIDDLE SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.




Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.BS 10.10	PC.BS 10.13	PC.BS 10.15
Yellow	PC.BS 20.10	PC.BS 20.13	PC.BS 20.15
Red	PC.BS 30.10	PC.BS 30.13	PC.BS 30.15
Blue	PC.BS 40.10	PC.BS 40.13	PC.BS 40.15
Green	PC.BS 50.10	PC.BS 50.13	PC.BS 50.15
Black	PC.BS 60.10	PC.BS 60.13	PC.BS 60.15
Orange	PC.BS 70.10	PC.BS 70.13	PC.BS 70.15
Purple	PC.BS 80.10	PC.BS 80.13	PC.BS 80.15
Bronze	PC.BS 90.10	PC.BS 90.13	PC.BS 90.15

# AIO SERIES - Valve Barrier Option

## 1/1 FULL SIZE



- AIO series - Valve Option containers are working with valve mechanism.
- AIO series containers filter systems can be detachable, valve parts can be changed easily.
- 9 Different colour options available.



Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.CS 10.10	PC.CS 10.13	PC.CS 10.15	PC.CS 10.20	PC.CS 10.26
Yellow	PC.CS 20.10	PC.CS 20.13	PC.CS 20.15	PC.CS 20.20	PC.CS 20.26
Red	PC.CS 30.10	PC.CS 30.13	PC.CS 30.15	PC.CS 30.20	PC.CS 30.26
Blue	PC.CS 40.10	PC.CS 40.13	PC.CS 40.15	PC.CS 40.20	PC.CS 40.26
Green	PC.CS 50.10	PC.CS 50.13	PC.CS 50.15	PC.CS 50.20	PC.CS 50.26
Black	PC.CS 60.10	PC.CS 60.13	PC.CS 60.15	PC.CS 60.20	PC.CS 60.26
Orange	PC.CS 70.10	PC.CS 70.13	PC.CS 70.15	PC.CS 70.20	PC.CS 70.26
Purple	PC.CS 80.10	PC.CS 80.13	PC.CS 80.15	PC.CS 80.20	PC.CS 80.26
Bronze	PC.CS 90.10	PC.CS 90.13	PC.CS 90.15	PC.CS 90.20	PC.CS 90.26




TW Surgical





# PLASMA MODEL STERILIZATION CONTAINERS





# PLASMA SERIES

## 1/2 HALF SIZE



- Plasma series containers can work only with the plasma sterilizer filters.



Colours

285 x 280 x 100mm

285 x 280 x 135mm

285 x 280 x 150mm

285 x 280 x 200mm

285 x 280 x 260mm

Gray

PC.AH 10.10

PC.AH 10.13

PC.AH 10.15


PC.AH 10.20

PC.AH 10.26

# PLASMA SERIES

## 3/4 MIDDLE SIZE

- Plasma series containers can work only with the plasma sterilizer filters.



Colours

465 x 280 x 100mm

465 x 280 x 135mm

465 x 280 x 150mm

Gray

PC.BH 10.10

PC.BH 10.13

PC.BH 10.15

# PLASMA SERIES

## 1/1 FULL SIZE



- Plasma series containers can work only with the plasma sterilizer filters.



Colours

580 x 280 x 100mm

580 x 280 x 135mm

580 x 280 x 150mm

580 x 280 x 200mm

580 x 280 x 260mm

Gray

PC.CH 10.10

PC.CH 10.13

PC.CH 10.15

PC.CH 10.20

PC.CH 10.26



**"Quality Assurance and an Innovative Approach"**



TW Surgical





TW Surgical

# REMOVAL MODEL STERILIZATION CONTAINERS




# REMOVAL SERIES

## 1/2 HALF SIZE



\*\*\*Attention!

- Removal containers can be use only for the transportation.



Colours	285 x 280 x 100mm	285 x 280 x 135mm	285 x 280 x 150mm	285 x 280 x 200mm	285 x 280 x 260mm
Gray	PC.AR 10.10	PC.AR 10.13	PC.AR 10.15	PC.AR 10.20	PC.AR 10.26
Yellow	PC.AR 20.10	PC.AR 20.13	PC.AR 20.15	PC.AS 20.20	PC.AR 20.26
Red	PC.AR 30.10	PC.AR 30.13	PC.AR 30.15	PC.AR 30.20	PC.AR 30.26
Blue	PC.AR 40.10	PC.AR 40.13	PC.AR 40.15	PC.AR 40.20	PC.AR 40.26
Green	PC.AR 50.10	PC.AR 50.13	PC.AR 50.15	PC.AR 50.20	PC.AR 50.26
Black	PC.AR 60.10	PC.AR 60.13	PC.AR 60.15	PC.AR 60.20	PC.AR 60.26
Orange	PC.AR 70.10	PC.AR 70.13	PC.AR 70.15	PC.AR 70.20	PC.AR 70.26
Purple	PC.AR 80.10	PC.AR 80.13	PC.AR 80.15	PC.AR 80.20	PC.AR 80.26
Bronze	PC.AR 90.10	PC.AR 90.13	PC.AR 90.15	PC.AR 90.20	PC.AR 90.26

# REMOVAL SERIES

## 3/4 MIDDLE SIZE

\*\*Attention!

- Removal containers can be use only for the transportation.



Colours	465 x 280 x 100mm	465 x 280 x 135mm	465 x 280 x 150mm
Gray	PC.BR 10.10	PC.BR 10.13	PC.BR 10.15
Yellow	PC.BR 20.10	PC.BR 20.13	PC.BR 20.15
Red	PC.BR 30.10	PC.BR 30.13	PC.BR 30.15
Blue	PC.BR 40.10	PC.BR 40.13	PC.BR 40.15
Green	PC.BR 50.10	PC.BR 50.13	PC.BR 50.15
Black	PC.BR 60.10	PC.BR 60.13	PC.BR 60.15
Orange	PC.BR 70.10	PC.BR 70.13	PC.BR 70.15
Purple	PC.BR 80.10	PC.BR 80.13	PC.BR 80.15
Bronze	PC.BR 90.10	PC.BR 90.13	PC.BR 90.15

# REMOVAL SERIES



## 1/1 FULL SIZE



TW Surgical

\*\*Attention!

- Removal containers can be use only for the transportation.




Colours	580 x 280 x 100mm	580 x 280 x 135mm	580 x 280 x 150mm	580 x 280 x 200mm	580 x 280 x 260mm
Gray	PC.CR 10.10	PC.CR 10.13	PC.CR 10.15	PC.CR 10.20	PC.CR 10.26
Yellow	PC.CR 20.10	PC.CR 20.13	PC.CR 20.15	PC.CR 20.20	PC.CR 20.26
Red	PC.CR 30.10	PC.CR 30.13	PC.CR 30.15	PC.CR 30.20	PC.CR 30.26
Blue	PC.CR 40.10	PC.CR 40.13	PC.CR 40.15	PC.CR 40.20	PC.CR 40.26
Green	PC.CR 50.10	PC.CR 50.13	PC.CR 50.15	PC.CR 50.20	PC.CR 50.26
Black	PC.CR 60.10	PCY.CR 60.13	PC.CR 60.15	PC.CR 60.20	PC.CR 60.26
Orange	PC.CR 70.10	PC.CR 70.13	PC.CR 70.15	PC.CR 70.20	PC.CR 70.26
Purple	PC.CR 80.10	PC.CR 80.13	PC.CR 80.15	PC.CR 80.20	PC.CR 80.26
Bronze	PC.CR 90.10	PC.CR 90.13	PC.CR 90.15	PC.CR 90.20	PC.CR 90.26

"Quality Assurance and an Innovative Approach"




TW Surgical





TW Surgical





TW Surgical

# DENTAL MODEL STERILIZATION CONTAINERS



# DENTAL SERIES STANDART



- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of Dental Size.
- 9 different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.



## BOTTOM NON-PERFORATED



Colours	310 x 190 x 30mm	310 x 190 x 50mm	310 x 190 x 70mm	310 x 190 x 100mm	310 x 190 x 130mm
Gray	PC.K 10.03	PC.K 10.05	PC.K 10.07	PC.K 10.10	PC.K 10.13
Yellow	PC.K 20.03	PC.K 20.05	PC.K 20.07	PC.K 20.10	PC.K 20.13
Red	PC.K 30.03	PC.K 30.05	PC.K 30.07	PC.K 30.10	PC.K 30.13
Blue	PC.K 40.03	PC.K 40.05	PC.K 40.07	PC.K 40.10	PC.K 40.13
Green	PC.K 50.03	PC.K 50.05	PC.K 50.07	PC.K 50.10	PC.K 50.13
Black	PC.K 60.03	PC.K 60.05	PC.K 60.07	PC.K 60.10	PC.K 60.13
Orange	PC.K 70.03	PC.K 70.05	PC.K 70.07	PC.K 70.10	PC.K 70.13
Purple	PC.K 80.03	PC.K 80.05	PC.K 80.07	PC.K 80.10	PC.K 80.13
Bronze	PC.K 90.03	PC.K 90.05	PC.K 90.07	PC.K 90.10	PC.K 90.13

# DENTAL SERIES STANDART



- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of Dental Size.
- 9 different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.



## BOTTOM PERFORATED




Colours	310 x 190 x 30mm	310 x 190 x 50mm	310 x 190 x 70mm	310 x 190 x 100mm	310 x 190 x 130mm
Gray	PC.L 10.03	PC.L 10.05	PC.L 10.07	PC.L 10.10	PC.L 10.13
Yellow	PC.L 20.03	PC.L 20.05	PC.L 20.07	PC.L 20.10	PC.L 20.13
Red	PC.L 30.03	PC.L 30.05	PC.L 30.07	PC.L 30.10	PC.L 30.13
Blue	PC.L 40.03	PC.L 40.05	PC.L 40.07	PC.L 40.10	PC.L 40.13
Green	PC.L 50.03	PC.L 50.05	PC.L 50.07	PC.L 50.10	PC.L 50.13
Black	PC.L 60.03	PC.L 60.05	PC.L 60.07	PC.L 60.10	PC.L 60.13
Orange	PC.L 70.03	PC.L 70.05	PC.L 70.07	PC.L 70.10	PC.L 70.13
Purple	PC.L 80.03	PC.L 80.05	PC.L 80.07	PC.L 80.10	PC.L 80.13
Bronze	PC.L 90.03	PC.L 90.05	PC.L 90.07	PC.L 90.10	PC.L 90.13

# DENTAL SERIES STANDART



TW Surgical

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of Dental Size.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.



## BOTTOM NON-PERFORATED



Colours	310 x 190 x 30mm	310 x 190 x 50mm	310 x 190 x 70mm	310 x 190 x 100mm	310 x 190 x 130mm
Gray	PC.K 00.03	PC.K 00.05	PC.K 00.07	PC.K 00.10	PC.K 00.13



## BOTTOM PERFORATED

Colours	310 x 190 x 30mm	310 x 190 x 50mm	310 x 190 x 70mm	310 x 190 x 100mm	310 x 190 x 130mm
Gray	PC.L 00.03	PC.L 00.05	PC.L 00.07	PC.L 00.10	PC.L 00.13

# DENTAL SERIES STANDART



## Lid

Colours 310 x 190mm

Gray PC.K 10.01

Yellow PC.K 20.01

Red PC.K 30.01

Blue PC.K 40.01

Green PC.K 50.01

Black PC.K 60.01


Orange PC.K 70.01

Purple PC.K 80.01

Bronze PC.K 90.01




TW Surgical






TW Surgical

# MINI MODEL STERILIZATION CONTAINERS



# MINI SERIES STANDART

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of Mini Size.
- 9 different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.




## BOTTOM NON-PERFORATED



Colours	300 x 140 x 30mm	300 x 140 x 50mm	300 x 140 x 70mm	300 x 140 x 100mm	300 x 140 x 130mm
<b>Gray</b>	PC.J 10.03	PC.J 10.05	PC.J 10.07	PC.J 10.10	PC.J 10.13
<b>Yellow</b>	PC.J 20.03	PC.J 20.05	PC.J 20.07	PC.J 20.10	PC.J 20.13
<b>Red</b>	PC.J 30.03	PC.J 30.05	PC.J 30.07	PC.J 30.10	PC.J 30.13
<b>Blue</b>	PC.J 40.03	PC.J 40.05	PC.J 40.07	PC.J 40.10	PC.J 40.13
<b>Green</b>	PC.J 50.03	PC.J 50.05	PC.J 50.07	PC.J 50.10	PC.J 50.13
<b>Black</b>	PC.J 60.03	PC.J 60.05	PC.J 60.07	PC.J 60.10	PC.J 60.13
<b>Orange</b>	PC.J 70.03	PC.J 70.05	PC.J 70.07	PC.J 70.10	PC.J 70.13
<b>Purple</b>	PC.J 80.03	PC.J 80.05	PC.J 80.07	PC.J 80.10	PC.J 80.13
<b>Bronze</b>	PC.J 90.03	PC.J 90.05	PC.J 90.07	PC.J 90.10	PC.J 90.13

# MINI SERIES STANDART

- Bottom Perforated and Bottom Non-Perforated options are available for Standart Series of Dental Size.
- 9 different colour options available.
- Single use paper filter, 100 cycles Textile filter and PTFE Filter options are available for Standart Seies.



## BOTTOM PERFORATED




Colours	300 x 140 x 30mm	300 x 140 x 50mm	300 x 140 x 70mm	300 x 140 x 100mm	300 x 140 x 130mm
Gray	PC.M 10.03	PC.M 10.05	PC.M 10.07	PC.M 10.10	PC.M 10.13
Yellow	PC.M 20.03	PC.M 20.05	PC.M 20.07	PC.M 20.10	PC.M 20.13
Red	PC.M 30.03	PC.M 30.05	PC.M 30.07	PC.M 30.10	PC.M 30.13
Blue	PC.M 40.03	PC.M 40.05	PC.M 40.07	PC.M 40.10	PC.M 40.13
Green	PC.M 50.03	PC.M 50.05	PC.M 50.07	PC.M 50.10	PC.M 50.13
Black	PC.M 60.03	PC.M 60.05	PC.M 60.07	PC.M 60.10	PC.M 60.13
Orange	PC.M 70.03	PC.M 70.05	PC.M 70.07	PC.M 70.10	PC.M 70.13
Purple	PC.M 80.03	PC.M 80.05	PC.M 80.07	PC.M 80.10	PC.M 80.13
Bronze	PC.M 90.03	PC.M 90.05	PC.M 90.07	PC.M 90.10	PC.M 90.13

# MINI SERIES STANDART



TW Surgical



## BOTTOM NON-PERFORATED

Colours	300 x 140 x 30mm	300 x 140 x 50mm	300 x 140 x 70mm	300 x 140 x 100mm	300 x 140 x 130mm
Gray	PC.J 00.03	PC.J 00.05	PC.J 00.07	PC.J 00.10	PC.J 00.13





## BOTTOM PERFORATED

Colours	300 x 140 x 30mm	300 x 140 x 50mm	300 x 140 x 70mm	300 x 140 x 100mm	300 x 140 x 130mm
Gray	PC.M 00.03	PC.M 00.05	PC.M 00.07	PC.M 00.10	PC.M 00.13

# MINI SERIES STANDART



TW Surgical



## Lid

Colours 300 x 140mm

Gray PC.J 10.01

Yellow PC.J 20.01

Red PC.J 30.01

Blue PC.J 40.01

Green PC.J 50.01

Black PC.J 60.01


Orange PC.J 70.01

Purple PC.J 80.01

Bronze PC.J 90.01



TW Surgical





TW Surgical

# PERFORATED SHEET BASKETS



# PERFORATED SHEET BASKET WITHOUT LID 1/2 HALF SIZE



Colours	255 x 245 x 30mm	255 x 245 x 50mm	255 x 245 x 70mm	255 x 245 x 100mm
Gray	PC.PA 10.03	PC.PA 10.05	PC.PA 10.07	PC.PA 10.10


# PERFORATED SHEET BASKET WITHOUT LID 3/4 MIDDLE SIZE




Colours	405 x 255 x 30mm	405 x 255 x 50mm	405 x 255 x 70mm	405 x 255 x 100mm
Gray	PC.PB 10.03	PC.PB 10.05	PC.PB 10.07	PC.PB 10.10

# PERFORATED SHEET BASKET WITHOUT LID

## 1/1 FULL SIZE



Colours	485 x 255 x 30mm	485 x 255 x 50mm	485 x 255 x 70mm	485 x 255 x 100mm
Gray	PC.PC 10.03	PC.PC 10.05	PC.PC 10.07	PC.PC 10.10



Colours	540 x 255 x 30mm	540 x 255 x 50mm	540 x 255 x 70mm	540 x 255 x 100mm
Gray	PC.PD 10.03	PC.PD 10.05	PC.PD 10.07	PC.PD 10.10

# PERFORATED SHEET BASKET WITH LID 1/2 HALF SIZE



TW Surgical



Colours	255 x 245 x 30mm	255 x 245 x 50mm	255 x 245 x 70mm	255 x 245 x 100mm
Gray	PC.KA 10.03	PC.KA 10.05	PC.KA 10.07	PC.KA 10.10

# PERFORATED SHEET BASKET WITH LID 3/4 MIDDLE SIZE




Colours	405 x 255 x 30mm	405 x 255 x 50mm	405 x 255 x 70mm	405 x 255 x 100mm
Gray	PC.KB 10.03	PC.KB 10.05	PC.KB 10.07	PC.KB 10.10

# PERFORATED SHEET BASKET WITH LID

## 1/1 FULL SIZE



Colours	485 x 255 x 30mm	485 x 255 x 50mm	485 x 255 x 70mm	485 x 255 x 100mm
Gray	PC.KC 10.03	PC.KC 10.05	PC.KC 10.07	PC.KC 10.10




Colours	540 x 255 x 30mm	540 x 255 x 50mm	540 x 255 x 70mm	540 x 255 x 100mm
Gray	PC.KD 10.03	PC.KD 10.05	PC.KD 10.07	PC.KD 10.10

# PERFORATED SHEET BASKET FOR DENTAL CONTAINER



TW Surgical




Colours	265 x 170 x 30mm	265 x 170 x 50mm	265 x 170 x 70mm	265 x 170 x 100mm
Gray	PC.DP 10.03	PC.DP 10.05	PC.DP 10.07	PC.DP 10.10

# PERFORATED SHEET BASKET FOR MINI CONTAINER




TW Surgical



Colours	265 x 125 x 30mm	265 x 125 x 50mm	265 x 125 x 70mm	265 x 125 x 100mm
Gray	PC.MP 10.03	PC.MP 10.05	PC.MP 10.07	PC.MP 10.10




TW Surgical





TW Surgical

# SILICONE MAT







TW Surgical

# SILICONE MAT 1/2 HALF SIZE

9 Colours Options




Colours	220 x 230 mm
Gray	PC.CM 10.13
Yellow	PC.CM 20.13
Red	PC.CM 30.13
Blue	PC.CM 40.13
Green	PC.CM 50.13
Black	PC.CM 60.13
Orange	PC.CM 70.13
Purple	PC.CM 80.13
Bronze	PC.CM 90.13




# SILICONE MAT 3/4 MIDDLE SIZE

9 Colours Options



Colours	380 x 230 mm
Gray	PC.CM 10.14
Yellow	PC.CM 20.14
Red	PC.CM 30.14
Blue	PC.CM 40.14
Green	PC.CM 50.14
Black	PC.CM 60.14
Orange	PC.CM 70.14
Purple	PC.CM 80.14
Bronze	PC.CM 90.14






# SILICONE MAT 1/1 FULL SIZE



TW Surgical

## 9 Colours Options



Colours	440 x 230 mm
Gray	PC.CM 10.15
Yellow	PC.CM 20.15
Red	PC.CM 30.15
Blue	PC.CM 40.15
Green	PC.CM 50.15
Black	PC.CM 60.15
Orange	PC.CM 70.15
Purple	PC.CM 80.15
Bronze	PC.CM 90.15


Colours	520 x 230 mm
Gray	PC.CM 10.16
Yellow	PC.CM 20.16
Red	PC.CM 30.16
Blue	PC.CM 40.16
Green	PC.CM 50.16
Black	PC.CM 60.16
Orange	PC.CM 70.16
Purple	PC.CM 80.16
Bronze	PC.CM 90.16

# SILICONE MATS FOR DENTAL


9 Colours Options



TW Surgical




Colours	250 x 150mm
Gray	PC.CM 10.10
Yellow	PC.CM 20.10
Red	PC.CM 30.10
Blue	PC.CM 40.10
Green	PC.CM 50.10
Black	PC.CM 60.10
Orange	PC.CM 70.10
Purple	PC.CM 80.10
Bronze	PC.CM 90.10




# SILICONE MATS FOR MINI

9 Colours Options



Colours	270 x 125mm
Gray	PC.CM 10.11
Yellow	PC.CM 20.11
Red	PC.CM 30.11
Blue	PC.CM 40.11
Green	PC.CM 50.11
Black	PC.CM 60.11
Orange	PC.CM 70.11
Purple	PC.CM 80.11
Bronze	PC.CM 90.11



# RACKS FOR ENDOSCOPY

\*\*\* Our accessories are custom designed. For more information feel free contact us.



TW Surgical

225 x 490 x 130 mm

PC.FS 10.01

Rack, for endoscopy

- Holds 10 endoscopic instruments
- Fits into a large size 20 cm deep container
- Wire basket large 50 mm
- Container large 200 mm
- Silicone mat 520 x 230 mm



225 x 490 x 130 mm

PC.FS 10.02

Rack, for endoscopy

- Holds 10 endoscopic instruments
- Fits into a large size 20 cm deep container
- Wire basket large 50 mm
- Container large 200 mm
- Silicone mat 520 x 230 mm



225 x 250 x 130 mm

PC.FS 10.03

Rack, for endoscopy

- Holds 10 endoscopic instruments
- Fits into a large size 20 cm deep container
- Wire basket large 50 mm
- Container large 200 mm
- Silicone mat 520 x 230 mm





Security Seal  
**1 Pack of 100 Pcs**  
 PC.KT 12.10



Basket Feet  
 PC.FT.10.0  
 1cm  $\varnothing$  2cm



Basket Feet  
 PC.FT.10.02  
 2cm  $\varnothing$  2cm

## PAPER LABEL



with indicator  
**1 Pack of 100 Pcs**  
 PC.ET 11.10



with indicator  
**1 Pack of 100 Pcs**  
 PC.ET 11.11


## ALUMINIUM LABEL



# ACCESSORIES




TW Surgical



Filter Retainer


PC.FT 10.10



Filter for Plasma Sterilization  
100 Cycles

PC.PF 10.03

**1 Pack of 100 Pcs**



Filter for Plasma Sterilization  
100 Cycles

PC.PF 10.04

**1 Pack of 100 Pcs**




Textile Filter  
100 Cycles

PC.BF 10.04



Textile Filter  
100 Cycles


PC.BF 10.03



Paper Filter  
Single Use


PC.KF 10.04

**1 Pack of 100 Pcs**



PTFE Filter  
2,000 Cycles

PC.TF 13.03



Paper Filter  
Single Use

PC.KF 10.03

**1 Pack of 100 Pcs**


# SILICONE FIXATION SYSTEM RAILS

\*\*\* For more information feel free contact us.




TW Surgical

with rail, screw and nut




- PC.FS 30.01
- 250 x 45 mm
- PC.FS 31.01
- 152 x 45 mm
- PC.FS 32.01
- 125 x 45 mm
- PC.FS 33.01
- 97 x 45 mm

with rail, screw and nut



- PC.FS 30.02
- 238 x 45 mm
- PC.FS 31.02
- 158 x 45 mm
- PC.FS 32.02
- 112 x 45 mm
- PC.FS 33.02
- 77 x 45 mm

with rail, screw and nut




- PC.FS 30.03
- 238 x 45 mm
- PC.FS 31.03
- 150 x 45 mm
- PC.FS 32.03
- 119 x 45 mm
- PC.FS 33.03
- 87 x 45 mm

with rail, screw and nut




- PC.FS 30.04
- 250 x 45 mm
- PC.FS 31.04
- 150 x 45 mm
- PC.FS 32.04
- 120 x 45 mm
- PC.FS 33.04
- 90 x 45 mm

Instruments holder, Ø 3 mm With rail, screw and nut




- PC.FS 30.05
- 178 x 24 mm

Instruments holder, Ø 3 to 8 mm With rail, screw and nut




- PC.FS 30.06
- 178 x 24 mm

Instruments holder, Ø 5 mm With rail, screw and nut




- PC.FS 30.07
- 73 x 24 mm

Instruments holder, Ø 3 mm With rail, screw and nut



- PC.FS 30.08
- 73 x 24 mm

with rail, screw and nut



- PC.FS 30.09
- 121 x 22 mm

# SILICONE FIXATION SYSTEM RAILS

\*\*\* For more information feel free contact us.




TW Surgical

Instruments holder, Ø 3 to 15 mm With rail, screw and nut




PC.FS 30.10  
178 x 34 mm

with rail, screw and nut



PC.FS 30.11  
123 x 22 mm

Instruments holder, Ø 5 mm With rail, screw and nut



PC.FS 30.12  
178 x 24 mm

Ø 8 mm With rail, screw and nut



PC.FS 30.13  
25.5 x 25 mm

Ø 10,5 mm With rail, screw and nut




PC.FS 30.14  
25.5 x 25 mm

Ø 13 mm With rail, screw and nut



PC.FS 30.15  
25.5 x 25 mm

Endoscope holder, double One side Ø 12 mm open One side from Ø 3 to Ø 15 mm closed With rail, screw and nut



PC.FS 30.16  
72 x 34 mm

Ø 5 mm With rail, screw and nut



PC.FS 30.17  
23.5 x 25 mm

Instruments holder from Ø 3 to Ø 15 mm With rail, screw and nut



PC.FS 30.18  
72 X 34 mm

Instruments holder, from Ø 3 to Ø 8 mm With rail, screw and nut




PC.FS 30.19  
72 x 24 mm

Endoscope holder, double one side from Ø 3 to Ø 8 mm With rail, screw and nut



PC.FS 30.20  
38 x 24 mm

Endoscope holder, double One side Ø 12 mm open One side from Ø 3 to Ø 8 mm closed With rail, screw and nut



PC.FS 30.21  
72 x 34 mm

# SILICONE FIXATION SYSTEM RAILS

\*\*\* For more information feel free contact us.



TW Surgical


Endoscope holder,  $\varnothing$  15 mm closed  
With rail, screw and nut



PC.FS 30.22

38 x 33 mm


Endoscope holder,  $\varnothing$  12 mm closed  
With rail, screw and nut



PC.FS 30.23

38 x 33 mm

Endoscope holder, single one side  $\varnothing$   
12 mm open one side from  $\varnothing$  3 to  $\varnothing$   
15 mm With rail, screw and nut



PC.FS 30.24

28 x 34 mm

Endoscope holder,  $\varnothing$  12 mm closed  
With rail, screw and nut



PC.FS 30.25

38 x 33 mm

Endoscope holder,  $\varnothing$  8 mm closed  
With rail, screw and nut



PC.FS 30.26

38 x 33 mm


Endoscope holder,  $\varnothing$  6 mm closed  
With rail, screw and nut



PC.FS 30.27

38 x 33 mm

Endoscope holder,  $\varnothing$  4 mm closed With  
rail, screw and nut



PC.FS 30.28

72 x 34 mm