



TW Surgical



**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com



Bottom non perforated



Lid colour

Dimensions

**600 x 285 x 110 mm**

silver	TWS 31379
yellow	TWS 31380Y
red	TWS 31381 R
blue	TWS 31382 BU
green	TWS 31383 G
black	TWS 31384 BA

110



**600 x 285 x 135 mm**

silver	TWS 31385
yellow	TWS 31386Y
red	TWS 31387 R
blue	TWS 31388 BU
green	TWS 31389 G
black	TWS 31390 BA

135



**600 x 285 x 150 mm**

silver	TWS 31391
yellow	TWS 31392Y
red	TWS 31393 R
blue	TWS 31394 BU
green	TWS 31395 G
black	TWS 31396 BA

150



**600 x 285 x 200 mm**

silver	TWS 31397
yellow	TWS 313980 Y
red	TWS 31399 R
blue	TWS 31400 BU
green	TWS 31401 G
black	TWS 31402 BA

200



**600 x 285 x 260 mm**

silver	TWS 31403
yellow	TWS 31404Y
red	TWS 31405 R
blue	TWS 31406 BU
green	TWS 31407 G
black	TWS 31408 BA

260



Container bottom




Bottom non perforated



Height	Dimensions
	<b>600 x 272 mm</b>
96 mm	<b>TWS 31409</b>
116 mm	<b>TWS 31410</b>
138 mm	<b>TWS 31412</b>
190 mm	<b>TWS 31413</b>
235 mm	<b>TWS 31414</b>

TAENIA Container lid




Lid colour	Dimensions
	<b>585 x 285 mm</b>
silver	<b>TWS 31415</b>
yellow	<b>TWS 31416</b>
red	<b>TWS 31417</b>
blue	<b>TWS 31418</b>
green	<b>TWS 31419</b>
black	<b>TWS 31420</b>

Wire basket




Height	Dimensions	Dimensions
	<b>480 x 245 mm</b>	<b>540 x 245 mm</b>
30 mm	<b>TWS 31421</b>	<b>TWS 31425</b>
50 mm	<b>TWS 31422</b>	<b>TWS 31426</b>
70 mm	<b>TWS 31423</b>	<b>TWS 31427</b>
100 mm	<b>TWS 31424</b>	<b>TWS 31428</b>

Lid for wire basket



Height	Dimensions	Dimensions
	<b>485 x 245 mm</b>	<b>540 x 245 mm</b>
5 mm	<b>TWS 31429</b>	<b>TWS 31430</b>

Stainless steel wire basket, with feet



Height	Dimensions	Dimensions
	<b>485 x 253 mm</b>	<b>540 x 253 mm</b>
40 mm	<b>TWS 31431</b>	<b>TWS 31435</b>
60 mm	<b>TWS 31432</b>	<b>TWS 31436</b>
80 mm	<b>TWS 31433</b>	<b>TWS 31437</b>
110 mm	<b>TWS 31434</b>	<b>TWS 31438</b>

Lid for stainless steel wire basket, with lock




Height	Dimensions	Dimensions
	<b>487 x 256 mm</b>	<b>542 x 256 mm</b>
4 mm	<b>TWS 31439</b>	<b>TWS 31440</b>

Silicone mat




Colour	Dimensions	Dimensions
<b>blue</b>	<b>440 x 230 mm</b>	<b>520 x 230 mm</b>
	<b>TWS 31441</b>	<b>TWS 31442</b>

Silicone mesh mat



Dimensions	Dimensions
<b>485 x 250 mm</b>	<b>534 x 247 mm</b>
<b>TWS 31443</b>	<b>TWS 31444</b>

Wire base



Dimensions
<b>545 x 255 x 30 mm</b>
<b>TWS 31445</b>

Bottom non perforated



Lid colour

Dimensions

**475 x 285 x 110 mm**

silver	TWS 31446
yellow	TWS 31447Y
red	TWS 31448 R
blue	TWS 31449 BU
green	TWS 31450 G
black	TWS 31451 BA

110



**475 x 285 x 135 mm**

silver	TWS 31452
yellow	TWS 31453 Y
red	TWS 31454 R
blue	TWS 31455 BU
green	TWS 31456 G
black	TWS 31457 BA

135



**475 x 285 x 150 mm**

silver	TWS 31458
yellow	TWS 31459 Y
red	TWS 31460 R
blue	TWS 31461 BU
green	TWS 31462 G
black	TWS 31463 BA

150



**475 x 285 x 200 mm**

silver	TWS 31464
yellow	TWS 31465 Y
red	TWS 31466 R
blue	TWS 31467 BU
green	TWS 31468 G
black	TWS 31469 BA

200



Container bottom



Bottom non perforated




Height

Dimensions

	<b>475 x 272 mm</b>
96 mm	<b>TWS 31470</b>
116 mm	<b>TWS 31471</b>
138 mm	<b>TWS 31472</b>
185 mm	<b>TWS 31473</b>

TAENIA Container lid



Lid colour

Dimensions

	<b>460 x 285 mm</b>
silver	<b>TWS 31474</b>
yellow	<b>TWS 31475</b>
red	<b>TWS 31476</b>
blue	<b>TWS 31477</b>
green	<b>TWS 31478</b>
black	<b>TWS 31479</b>

Wire basket



Height

Dimensions

	<b>405 x 245 mm</b>
30 mm	<b>TWS 31480</b>
50 mm	<b>TWS 31481</b>
70 mm	<b>TWS 31482</b>
100 mm	<b>TWS 31483</b>

Lid for wire basket




Height

Dimensions

	<b>405 x 245 mm</b>
5 mm	<b>TWS 31484</b>

## Accessories for TAENIA Container Size 3/4

Stainless steel wire basket, with feet




Height	Dimensions
	<b>405 x 253 mm</b>
40 mm	<b>TWS 31485</b>
60 mm	<b>TWS 31486</b>
80 mm	<b>TWS 31487</b>
110 mm	<b>TWS 31488</b>

Lid for stainless steel wire basket, with lock




Height	Dimensions
	<b>407 x 256 mm</b>
4 mm	<b>TWS 31489</b>

Silicone mat



Colour	Dimensions
	<b>380 x 230 mm</b>
blue	<b>TWS 31490</b>

Silicone mesh mat




Dimensions
<b>395 x 245 mm</b>
<b>TWS 31491</b>

Wire base



Dimensions
<b>405 x 255 x 30 mm</b>
<b>TWS 31492</b>



Bottom non perforated



Lid colour

Dimensions

**310 x 280 x 110 mm**

silver	TWS 31493
yellow	TWS 31494Y
red	TWS 31495 R
blue	TWS 31496 BU
green	TWS 31497 G
black	TWS 31498 BA

110



**310 x 280 x 135 mm**

silver	TWS 31499
yellow	TWS 31500Y
red	TWS 31501 R
blue	TWS 31502 BU
green	TWS 31503 G
black	TWS 31504 BA

135



**310 x 280 x 150 mm**

silver	TWS 31505
yellow	TWS 31506Y
red	TWS 31507 R
blue	TWS 31508 BU
green	TWS 31509 G
black	TWS 31510 BA

150



**310 x 280 x 200 mm**

silver	TWS 31511
yellow	TWS 31512Y
red	TWS 31513 R
blue	TWS 31514 BU
green	TWS 31515 G
black	TWS 31516 BA

200



**310 x 280 x 260 mm**

silver	TWS 31517
yellow	TWS 31518Y
red	TWS 31519 R
blue	TWS 31520 BU
green	TWS 31521 G
black	TWS 31522 BA

260



Container bottom



Bottom non perforated



Height

Dimensions

**310 x 272 mm**

96 mm

**TWS 31523**

116 mm

**TWS 31524**

138 mm

**TWS 31525**


185 mm

**TWS 31526**

235 mm

**TWS 31527**

TAENIA Container lid



Lid colour

Dimensions

**295 x 280 mm**

silver

**TWS 31528**

yellow

**TWS 31529**

red

**TWS 31530**

blue

**TWS 31531**

green

**TWS 31532**

black

**TWS 31533**

Wire basket



Height

Dimensions

**246 x 246 mm**

30 mm

**TWS 31534**

50 mm

**TWS 31535**

70 mm

**TWS 31536**

100 mm

**TWS 31537**

Lid for wire basket



Height


Dimensions

**246 x 246 mm**

5 mm


**TWS 31538**

Stainless steel wire basket, with feet




Height	Dimensions
	<b>253 x 244 mm</b>
40 mm	<b>TWS 31539</b>
60 mm	<b>TWS 31540</b>
80 mm	<b>TWS 31541</b>
110 mm	<b>TWS 31542</b>

Lid for stainless steel wire basket, with lock



Height	Dimensions
	<b>256 x 247 mm</b>
4 mm	<b>TWS 31543</b>

Silicone mat




Colour	Dimensions
	<b>230 x 220 mm</b>
<b>blue</b>	<b>TWS 31544</b>

Silicone mesh mat



Dimensions
<b>250 x 240 mm</b>
<b>TWS 31545</b>

Wire base



Dimensions
<b>255 x 255 x 30 mm</b>
<b>TWS 31546</b>

Standard Container Size 1/1



Bottom non perforated

Bottom perforated



Lid colour

Dimensions

Dimensions

600 x 285 x 110 mm

600 x 285 x 110 mm

silver
yellow
red
blue
green
black

TWS 31547
TWS 31549Y
TWS 31550 R
TWS 31551 BU
TWS 31552 G
TWS 31553 BA

TWS 31548
TWS 31554Y
TWS 31555 R
TWS 31556 BU
TWS 31557 G
TWS 31558 BA

110



600 x 285 x 135 mm

600 x 285 x 135 mm

silver
yellow
red
blue
green
black

TWS 31559
TWS 31561Y
TWS 31562 R
TWS 31563 BU
TWS 31564 G
TWS 31565 BA

TWS 31560
TWS 31566 Y
TWS 31567 R
TWS 31568 BU
TWS 31569 G
TWS 31570 BA

135



600 x 285 x 150 mm

600 x 285 x 150 mm

silver
yellow
red
blue
green
black

TWS 31571
TWS 31572Y
TWS 31573 R
TWS 31574 BU
TWS 31575 G
TWS 31576 BA

TWS 31577
TWS 31578 Y
TWS 31579 R
TWS 31580 BU
TWS 31581 G
TWS 31582 BA

150



600 x 285 x 200 mm

600 x 285 x 200 mm

silver
yellow
red
blue
green
black

TWS 31584
TWS 31585Y
TWS 31586 R
TWS 31587 BU
TWS 31588 G
TWS 31589 BA

TWS 31583
TWS 31590 Y
TWS 31591 R
TWS 31592 BU
TWS 31593 G
TWS 31594 BA

200



600 x 285 x 260 mm

600 x 285 x 260 mm




silver
yellow
red
blue
green
black


TWS 31595
TWS 31602Y
TWS 31603 R
TWS 31604 BU
TWS 31605 G
TWS 31606 BA


TWS 31596
TWS 31597 Y
TWS 31598 R
TWS 31599 BU
TWS 31600 G
TWS 31601 BA


260




Container bottom		Bottom non perforated		Bottom perforated
		Height	Dimensions	Dimensions
				
			<b>600 x 272 mm</b>	<b>600 x 272 mm</b>
		96 mm	<b>TWS 31607</b>	<b>TWS 31615</b>
		116 mm	<b>TWS 31608</b>	<b>TWS 31616</b>
		138 mm	<b>TWS 31609</b>	<b>TWS 31612</b>
		190 mm	<b>TWS 31610</b>	<b>TWS 31613</b>
		235 mm	<b>TWS 31611</b>	<b>TWS 31614</b>


Container lid with filter holding device		Lid colour	Dimensions
			<b>585 x 285 mm</b>
	silver	<b>TWS 31617</b>	
	yellow	<b>TWS 31618</b>	
	red	<b>TWS 31619</b>	
	blue	<b>TWS 31620</b>	
	green	<b>TWS 31621</b>	
	black	<b>TWS 31622</b>	

Safety lid		Dimensions	
		Lid colour	Dimensions
			<b>606 x 288 mm</b>
		silver	<b>TWS 31623</b>
		yellow	<b>TWS 31624</b>
		red	<b>TWS 31625</b>
		blue	<b>TWS 31626</b>
		green	<b>TWS 31627</b>
		black	<b>TWS 31628</b>

Wire basket		Height	Dimensions	Dimensions
				<b>480 x 245 mm</b>
		30 mm	<b>TWS 31629</b>	<b>TWS 31633</b>
		50 mm	<b>TWS 31630</b>	<b>TWS 31634</b>
		70 mm	<b>TWS 31631</b>	<b>TWS 31635</b>
		100 mm	<b>TWS 31632</b>	<b>TWS 31636</b>

Lid for wire basket		Height	Dimensions	Dimensions
				<b>485 x 245 mm</b>
		5 mm	<b>TWS 31637</b>	<b>TWS 31638</b>

Stainless steel wire basket, with feet



Height

Dimensions

Dimensions

40 mm

**485 x 253 mm**

**540 x 253 mm**

60 mm

**TWS 31639**

**TWS 31643**

80 mm

**TWS 31640**

**TWS 31645**

110 mm

**TWS 31641**

**TWS 31647**

**TWS 31642**

**TWS 31648**

Lid for stainless steel wire basket, with lock



Height

Dimensions

Dimensions

4 mm

**487 x 256 mm**

**542 x 256 mm**

**TWS 31649**

**TWS 31650**

Silicone mat



Colour

Dimensions

Dimensions

blue


**440 x 230 mm**

**520 x 230 mm**

**TWS 31651**

**TWS 31652**

Silicone mesh mat



Dimensions

Dimensions

**485 x 250 mm**

**534 x 247 mm**

**TWS 31653**

**TWS 31654**

Wire base



Dimensions

**545 x 255 x 30 mm**

**TWS 31656**

Standard Container Size 3/4



Bottom non perforated

Bottom perforated



Lid colour

Dimensions

Dimensions

475 x 285 x 110 mm

475 x 285 x 110 mm

silver
yellow
red
blue
green
black

TWS 31657
TWS 31658Y
TWS 31659 R
TWS 31660 BU
TWS 31661 G
TWS 31662 BA

TWS 31663
TWS 31664Y
TWS 31665 R
TWS 31666 BU
TWS 31667 G
TWS 31668 BA

475 x 285 x 135 mm

475 x 285 x 135 mm

silver
yellow
red
blue
green
black

TWS 31669
TWS 31672 Y
TWS 31673 R
TWS 31674 BU
TWS 31675 G
TWS 31676 BA

TWS 31670
TWS 31677Y
TWS 31678 R
TWS 31679 BU
TWS 31680 G
TWS 31682 BA

475 x 285 x 150 mm

475 x 285 x 150 mm

silver
yellow
red
blue
green
black

TWS 31684
TWS 31685 Y
TWS 31686 R
TWS 31687 BU
TWS 31688 G
TWS 31689 BA

TWS 31683
TWS 31690Y
TWS 31691 R
TWS 31692 BU
TWS 31693 G
TWS 31694 BA

475 x 285 x 200 mm

475 x 285 x 200 mm

silver
yellow
red
blue
green
black

TWS 31695
TWS 31697Y
TWS 31706 R
TWS 31705 BU
TWS 31704 G
TWS 31703 BA

TWS 31696
TWS 31698Y
TWS 31699 R
TWS 31700 BU
TWS 31701 G
TWS 31702 BA

110



135






150




200



Container bottom		Bottom non perforated		Bottom perforated
		Height	Dimensions	Dimensions
				
			<b>475 x 272 mm</b>	<b>475 x 272 mm</b>
		96 mm	<b>TWS 31707</b>	<b>TWS 31714</b>
		116 mm	<b>TWS 31708</b>	<b>TWS 31711</b>
		138 mm	<b>TWS 31709</b>	<b>TWS 31712</b>
		185 mm	<b>TWS 31710</b>	<b>TWS 31713</b>


Container lid with filter holding device		Lid colour	Dimensions
			<b>460 x 285 mm</b>
	silver	<b>TWS 31715</b>	
	yellow	<b>TWS 31716</b>	
	red	<b>TWS 31717</b>	
	blue	<b>TWS 31718</b>	
	green	<b>TWS 31719</b>	
	black	<b>TWS 31720</b>	

Safety lid		Dimensions	
		<b>483 x 290 mm</b>	
	silver	<b>TWS 31721</b>	
	yellow	<b>TWS 31726</b>	
	red	<b>TWS 31722</b>	
	blue	<b>TWS 31723</b>	
	green	<b>TWS 31724</b>	
	black	<b>TWS 31725</b>	

Wire basket		Height	Dimensions
			<b>405 x 245 mm</b>
	30 mm	<b>TWS 31727</b>	
	50 mm	<b>TWS 31728</b>	
	70 mm	<b>TWS 31729</b>	
	100 mm	<b>TWS 31730</b>	


Lid for wire basket		Height	Dimensions
			<b>405 x 245 mm</b>
	5 mm	<b>TWS 31731</b>	

Stainless steel wire basket, with feet




Height	Dimensions
	<b>405 x 253 mm</b>
40 mm	<b>TWS 31732</b>
60 mm	<b>TWS 31733</b>
80 mm	<b>TWS 31734</b>
110 mm	<b>TWS 31735</b>

Lid for stainless steel wire basket, with lock




Height	Dimensions
	<b>407 x 256 mm</b>
4 mm	<b>TWS 31736</b>

Silicone mat




Colour	Dimensions
	<b>380 x 230 mm</b>
blue	<b>TWS 31737</b>

Silicone mesh mat



Dimensions
<b>395 x 245 mm</b>
<b>TWS 31738</b>

Wire base



Dimensions
<b>405 x 255 x 30 mm</b>
<b>TWS 31739</b>

Standard Container Size 1/2



Bottom non perforated



Bottom perforated



110



Lid colour

Dimensions

Dimensions

	310 x 280 x 110 mm	310 x 280 x 110 mm
silver	TWS 31740	TWS 31746
yellow	TWS 31741 Y	TWS 31747 Y
red	TWS 31742 R	TWS 31748 R
blue	TWS 31743 BU	TWS 31749 BU
green	TWS 31744 G	TWS 31750 G
black	TWS 31745 BA	TWS 31751 BA

135



	310 x 280 x 135 mm	310 x 280 x 135 mm
silver	TWS 31752	TWS 31758
yellow	TWS 31753 Y	TWS 31759 Y
red	TWS 31754 R	TWS 31760 R
blue	TWS 31755 BU	TWS 31761 BU
green	TWS 31756 G	TWS 31762 G
black	TWS 31757 BA	TWS 31763 BA

150



	310 x 280 x 150 mm	310 x 280 x 150 mm
silver	TWS 31764	TWS 31770
yellow	TWS 31765 Y	TWS 31771 Y
red	TWS 31766 R	TWS 31772 R
blue	TWS 31767 BU	TWS 31773 BU
green	TWS 31768 G	TWS 31774 G
black	TWS 31769 BA	TWS 31775 BA

200






	310 x 280 x 200 mm	310 x 280 x 200 mm
silver	TWS 31776	TWS 31777
yellow	TWS 31778 Y	TWS 31783 Y
red	TWS 31779 R	TWS 31784 R
blue	TWS 31780 BU	TWS 31785 BU
green	TWS 31781 G	TWS 31786 G
black	TWS 31782 BA	TWS 31787 BA

260



	310 x 280 x 260 mm	310 x 280 x 260 mm
silver	TWS 31788	TWS 31789
yellow	TWS 31795 Y	TWS 31790 Y
red	TWS 31796 R	TWS 31791 R
blue	TWS 31798 BU	TWS 31792 BU
green	TWS 31797 G	TWS 31793 G
black	TWS 31799 BA	TWS 31794 BA

Container bottom		Height	Bottom non perforated	Bottom perforated
			Dimensions	Dimensions
				
			<b>310 x 272 mm</b>	<b>310 x 272 mm</b>
		96 mm	<b>TWS 31800</b>	<b>TWS 31805</b>
		116 mm	<b>TWS 31801</b>	<b>TWS 31806</b>
		138 mm	<b>TWS 31802</b>	<b>TWS 31807</b>
		185 mm	<b>TWS 31803</b>	<b>TWS 31808</b>
		235 mm	<b>TWS 31804</b>	<b>TWS 31809</b>


Container lid with filter holding device		Lid colour	Dimensions
			<b>295 x 280 mm</b>
		silver	<b>TWS 31810</b>
		yellow	<b>TWS 31811</b>
		red	<b>TWS 31812</b>
		blue	<b>TWS 31813</b>
		green	<b>TWS 31814</b>
		black	<b>TWS 31815</b>

Safety lid			Dimensions
			<b>317 x 285 mm</b>
		silver	<b>TWS 31816</b>
		yellow	<b>TWS 31817</b>
		red	<b>TWS 31818</b>
		blue	<b>TWS 31819</b>
		green	<b>TWS 31820</b>
		black	<b>TWS 31821</b>

Wire basket		Height	Dimensions
			<b>246 x 246 mm</b>
		30 mm	<b>TWS 31822</b>
		50 mm	<b>TWS 31823</b>
		70 mm	<b>TWS 31824</b>
		100 mm	<b>TWS 31825</b>

Lid for wire basket		Height	Dimensions
			<b>246 x 246 mm</b>
		5 mm	<b>TWS 31826</b>

Stainless steel wire basket, with feet




Height

Dimensions

	<b>253 x 244 mm</b>
40 mm	<b>TWS 31829</b>
60 mm	<b>TWS 31830</b>
80 mm	<b>TWS 31827</b>
110 mm	<b>TWS 31828</b>

Lid for stainless steel wire basket, with lock




Height

Dimensions

	<b>256 x 247 mm</b>
4 mm	<b>TWS 31831</b>

Silicone mat



Colour

Dimensions

	<b>230 x 220 mm</b>
blue	<b>TWS 31832</b>


Silicone mesh mat



Dimensions

<b>250 x 240 mm</b>
<b>TWS 31833</b>

Wire base



Dimensions

<b>255 x 255 x 30 mm</b>
<b>TWS 31834</b>

# Standard Container Flat



TW Surgical



Bottom non perforated



Bottom perforated



Lid colour

Dimensions

Dimensions

**295 x 280 x 75 mm**

**295 x 280 x 75 mm**

silver	<b>TWS 31835</b>	<b>TWS 31836</b>
yellow	<b>TWS 31837Y</b>	<b>TWS 31842Y</b>
red	<b>TWS 31838 R</b>	<b>TWS 31843 R</b>
blue	<b>TWS 31839 BU</b>	<b>TWS 31844BU</b>
green	<b>TWS 31840 G</b>	<b>TWS 31845 G</b>
black	<b>TWS 31841 BA</b>	<b>TWS 31846 BA</b>

**295 x 280 x 90 mm**

**295 x 280 x 90 mm**

silver	<b>TWS 31847</b>	<b>TWS 31853</b>
yellow	<b>TWS 31848Y</b>	<b>TWS 31854Y</b>
red	<b>TWS 31849 R</b>	<b>TWS 31855 R</b>
blue	<b>TWS 31850 BU</b>	<b>TWS 31856 BU</b>
green	<b>TWS 31851 G</b>	<b>TWS 31857 G</b>
black	<b>TWS 31852 BA</b>	<b>TWS 31858 BA</b>

75







90





Treasure Wright (U.K.) Ltd.


128-160 City Road, London, EC1V 2NX, United Kingdom  
info@twrightmedical.com

Container bottom		Bottom non perforated		Bottom perforated
		Height	Dimensions	Dimensions
				
			<b>295 x 270 mm</b>	<b>295 x 270 mm</b>
		60 mm	<b>TWS 31859</b>	<b>TWS 31861</b>
		72 mm	<b>TWS 31860</b>	<b>TWS 31862</b>


Container lid with filter holding device		Lid colour	Dimensions
			<b>295 x 280 mm</b>
	silver	<b>TWS 31863</b>	
	yellow	<b>TWS 31866</b>	
	red	<b>TWS 31867</b>	
	blue	<b>TWS 31868</b>	
	green	<b>TWS 31864</b>	
	black	<b>TWS 31865</b>	

Safety lid		317 x 285 mm	
		silver	<b>TWS 31869</b>
	yellow	<b>TWS 31870</b>	
	red	<b>TWS 31872</b>	
	blue	<b>TWS 31873</b>	
	green	<b>TWS 31874</b>	
	black	<b>TWS 31875</b>	

Wire basket		Height	Dimensions
			<b>246 x 246 mm</b>
	30 mm	<b>TWS 31871</b>	
	50 mm	<b>TWS 31876</b>	

Lid for wire basket		Height	Dimensions
			<b>246 x 246 mm</b>
	5 mm	<b>TWS 31877</b>	

Stainless steel wire basket, with feet



Height

Dimensions

40 mm


**253 x 244 mm**

60 mm

**TWS 31878**

**TWS 31879**

Lid for stainless steel wire basket, with lock



Height


Dimensions

4 mm

**256 x 247 mm**

**TWS 31880**

Silicone mat



Colour

Dimensions

blue

**230 x 220 mm**

**TWS 31881**

Silicone mesh mat




Dimensions

**250 x 240 mm**

**TWS 31882**

Wire base



Dimensions

**255 x 255 x 30 mm**

**TWS 31883**

# Standard Container Mini



Bottom non perforated

Bottom perforated



Lid colour

Dimensions

Dimensions

**315 x 140 x 45 mm**

**315 x 140 x 45 mm**

silver	TWS 31884	TWS 31890
yellow	TWS 31885Y	TWS 31891Y
red	TWS 31886 R	TWS 31892 R
blue	TWS 31887 BU	TWS 31893 BU
green	TWS 31888 G	TWS 31894 G
black	TWS 31889 BA	TWS 31895 BA

45



**315 x 140 x 70 mm**

**315 x 140 x 70 mm**

silver	TWS 31896	TWS 31902
yellow	TWS 31897Y	TWS 31903 Y
red	TWS 31898 R	TWS 31904 R
blue	TWS 31899 BU	TWS 31905 BU
green	TWS 31900 G	TWS 31906 G
black	TWS 31901 BA	TWS 31907 BA

70



**315 x 140 x 100 mm**

**315 x 140 x 100 mm**



silver	TWS 31908	TWS 31914
yellow	TWS 31909Y	TWS 31915Y
red	TWS 31910 R	TWS 31916 R
blue	TWS 31911 BU	TWS 31917 BU
green	TWS 31912 G	TWS 31918 G
black	TWS 31913 BA	TWS 31919 BA

100



Container bottom



	Bottom non perforated	Bottom perforated
		
Height		
	Dimensions	Dimensions
	<b>315 x 132 mm</b>	<b>315 x 132 mm</b>
32 mm	<b>TWS 31920</b>	<b>TWS 31923</b>
63 mm	<b>TWS 31921</b>	<b>TWS 31924</b>
92 mm	<b>TWS 31922</b>	<b>TWS 31925</b>

Container lid with filter holding device



Lid colour	Dimensions
	<b>315 x 140 mm</b>
silver	<b>TWS 31926</b>
yellow	<b>TWS 31927</b>
red	<b>TWS 31928</b>
blue	<b>TWS 31929</b>
green	<b>TWS 31930</b>
black	<b>TWS 31931</b>

Wire basket




Height	Dimensions
	<b>268 x 117 mm</b>
43 mm	<b>TWS 31932</b>

Wire basket with Lid



Height	Dimensions
	<b>268 x 127 mm</b>
47 mm	<b>TWS 31933</b>

Stainless steel wire basket, with feet



Height


55 mm

Dimensions

**265 x 125 mm**

**TWS 31934**

Foot for Mini and Special stainless steel trays, height 5 mm



Item number

**TWS 31935**

Lid for stainless steel wire basket, with lock



Dimensions

**267 x 127 mm**

**TWS 31936**

Silicone mat



Colour

blue

Dimensions

**270 x 125 mm**

**TWS 31937**



Bottom non perforated

Bottom perforated



55



Lid colour	Dimensions	Dimensions
	<b>325 x 192 x 55 mm</b>	<b>325 x 192 x 55 mm</b>
silver	TWS 31938	TWS 31945
yellow	TWS 31940 Y	TWS 31946 Y
red	TWS 31941 R	TWS 31947 R
blue	TWS 31942 BU	TWS 31948 BU
green	TWS 31943 G	TWS 31949 G
black	TWS 31944 BA	TWS 31950 BA

70



Lid colour	Dimensions	Dimensions
	<b>325 x 192 x 70 mm</b>	<b>325 x 192 x 70 mm</b>
silver	TWS 31951	TWS 31957
yellow	TWS 31952 Y	TWS 31958 Y
red	TWS 31953 R	TWS 31959 R
blue	TWS 31954 BU	TWS 31960 BU
green	TWS 31955 G	TWS 31961 G
black	TWS 31956 BA	TWS 31962 BA

100



Lid colour	Dimensions	Dimensions
	<b>325 x 192 x 100 mm</b>	<b>325 x 192 x 100 mm</b>
silver	TWS 31963	TWS 31969
yellow	TWS 31964 Y	TWS 31970 Y
red	TWS 31965 R	TWS 31971 R
blue	TWS 31966 BU	TWS 31972 BU
green	TWS 31967 G	TWS 31973 G
black	TWS 31968 BA	TWS 31974 BA

135



Lid colour	Dimensions	Dimensions
	<b>325 x 192 x 135 mm</b>	<b>325 x 192 x 135 mm</b>
silver	TWS 31975	TWS 31976
yellow	TWS 31977 Y	TWS 31982 Y
red	TWS 31978 R	TWS 31983 R
blue	TWS 31979 BU	TWS 31984 BU
green	TWS 31980 G	TWS 31985 G
black	TWS 31981 BA	TWS 31986 BA

Container bottom



Height

43 mm
59 mm
102 mm
122 mm

Bottom non perforated



Dimensions

**320 x 186 mm**

**TWS 31987**

**TWS 31988**

**TWS 31989**

**TWS 31990**

Bottom perforated



Dimensions

**320 x 186 mm**


**TWS 31991**

**TWS 31992**

**TWS 31993**

**TWS 31994**

Container lid with filter holding device



Lid colour

silver
yellow
red
blue
green
black

Dimensions

**325 x 192 mm**

**TWS 31995**

**TWS 31996**


**TWS 31997**

**TWS 31998**

**TWS 31999**

**TWS 32000**

Stainless steel wire basket, with feet



Height

45 mm
105 mm


Dimensions

**275 x 175 mm**

**TWS 32001**

**TWS 32002**

Lid for stainless steel wire basket, with lock



Height

4 mm
------

Dimensions

**277 x 178 mm**

**TWS 32003**

Foot for Mini and Special stainless steel trays, height 5 mm

Item number

**TWS 32004**



Silicone mat


Colour

Dimensions

**260 x 180 mm**

blue

**TWS 32005**



Endoscopy baskets with silicone holders for endoscopes




Dimensions	1 Endoscope	2 Endoscopes
290 x 80 x 55 mm	<b>TWS 32006</b>	<b>TWS 32007</b>
460 x 80 x 55 mm	<b>TWS 32008</b>	<b>TWS 32009</b>
670 x 80 x 55 mm	<b>TWS 32010</b>	

Endoscopy baskets with silicone holders for a variable numbers of endoscopes



Dimensions	
290 x 160 x 52 mm	<b>TWS 32011</b>
310 x 150 x 77 mm	<b>TWS 32012</b>
460 x 160 x 52 mm	<b>TWS 32013</b>
460 x 225 x 52 mm	<b>TWS 32014</b>

Sterilization baskets for a variable numbers of endoscopes / URS with adjustable holding device




Dimensions	
580 x 160 x 52 mm	<b>TWS 32015</b>
640 x 150 x 77 mm	<b>TWS 32016</b>

! Delivery without endoscope

## Endoscopy Baskets, Scope Containers

laparoscope and light cable, complete with racks and lid



Dimensions

**506 x 172 x 56 mm**

**TWS 32017 NIR**

Basket for Cystoscope/Hysteroscope and light cable, complete with racks and lid



Dimensions

**406 x 172 x 56 mm**

**TWS 32018**

Alu scope container, without filters, only for transport and storage, with endoscope holding devices



Dimensions

275 x 70 x 65 mm


**TWS 32019**

470 x 82 x 74 mm

**TWS 32020**

**Treasure Wright (U.K.) Ltd.**


Sterilization basket for 2 endoscopes and light cable



Dimensions

**492 x 134 x 65 mm**

**TWS 32021**



! Delivery without endoscope and light cable

Sterilizing container Endo, complete



Bottom non perforated



Bottom perforated



Lid colour

silver

Dimensions

**532 x 160 x 100 mm**

**TWS 32022**

Dimensions

**532 x 160 x 100 mm**

**TWS 32023**

Container bottom



Height

88 mm

Dimensions

**528 x 145 mm**

**TWS 32024**

Dimensions

**528 x 145 mm**

**TWS 32025**

Container lid with filter holding device



Lid colour


silver

Dimensions

**532 x 160 mm**

**TWS 32026**

Baskets for Laparoscopy Basic Set, complete with racks and lid (consisting of Tray 01 + 02)




Dimensions

**540 x 243 x 116 mm**

**TWS 32027**

Basket, complete with racks and lid



Dimensions

**540 x 243 x 66 mm**

**TWS 32028**

for Resectoscopy/Cystoscopy Set


Dimensions

**540 x 243 x 66 mm**

**TWS 32029**

for Hysteroscopy Set

Endoscopy rack with silicone holders for 10 instruments Ø 5.0 | 10.0 mm




Dimensions

**485 x 235 mm**

**TWS 32030**

Endoscopy wire basket, complete with 2 racks and lid




! Delivery without instruments

Dimensions

**530 x 250 x 200 mm**

**TWS 32031**

Insert for Resectoscopy Set



Item number

**TWS 32032**

TAENIA BioStop barrier disc, PPSU



Item number

**TWS 32033**

Disposable paper filter for standard container size 1/1, 3/4, 1/2 und Flat




Dimensions

**Ø 184 mm**

**TWS 32034**

Disposable polypropylen filter for standard container size 1/1, 3/4, 1/2 und Flat, for plasma sterilization



Dimensions

**Ø 184 mm**

**TWS 32035**

PTFE Permanent filter for standard container size 1/1, 3/4, 1/2 und Flat



max. Cycles

Dimensions

**Ø 190 mm**

1,200

**TWS 32036**

2,000

**TWS 32037**

Filter holder for Standard container size 1/1, 3/4, 1/2 und Flat



Dimensions

**Ø 186 mm**

**TWS 32038**

Dispenser for disposable paper filters



Dimensions

**200 x 200 x 64 mm**

**TWS 32039**

Disposable paper filter for container Mini




Dimensions

**235 x 120 mm**

**TWS 32040**

Disposable polypropylen filter for container Mini, for plasma sterilization




Dimensions

**235 x 120 mm**

**TWS 32041**

PTFE Permanent filter for container Mini



Cycles

Dimensions

**215 x 100 mm**

1,200

**TWS 32042**

2,000

**TWS 32043**

Disposable paper filter for container Special




Dimensions

**235 x 165 mm**

**TWS 32044**

Disposable polypropylen filter for container Special, for plasma sterilization




Dimensions

**235 x 165 mm**

**TWS 32045**

PTFE Permanent filter for container Special



max. Cycles    Dimensions

**215 x 145 mm**

**2,000**

**TWS 32046**

Disposable paper filter for Endo container




Dimensions

**235 x 105 mm**

**TWS 32047**

Disposable polypropylen filter for Endo container, for plasma sterilization



Dimensions

**235 x 105 mm**

**TWS 32048**


PTFE Permanent filter for Endo Container



max. Cycles    Dimensions

	<b>215 x 95 mm</b>
<b>2,000</b>	<b>TWS 32049</b>

Indicator label for container size 1/1, 3/4 and 1/2



Item number

**TWS 32050**



Indicator label for container Mini, Special and Flat



Item number


**TWS 32051**

Indicator label for container size 1/1, 3/4 and 1/2, for plasma sterilization




Item number

**TWS 32052**




Indicator label for container Mini, Special and Flat, for plasma sterilization




Item number

**TWS 32053**



Auto lock mechanism for container size 1/1, 3/4 and 1/2



Item number

**TWS 32054**



Order quantity: 2 pieces

Security seal without indicator



Colour

Item number

blue

**TWS 32055**

red

**TWS 32056**



Order quantity: 100 pieces

Security seal with indicator



Colour

Item number

blue

**TWS 32057**

red

**TWS 32058**



Order quantity: 100 pieces

Security seal with indicator, for plasma sterilization



Colour

Item number

blue

**TWS 32059**



Order quantity: 100 pieces

Identification labels



Colour	Item number
silver	<b>TWS 32060</b>
yellow	<b>TWS 32062</b>
red	<b>TWS 32063</b>
blue	<b>TWS 32064</b>
green	<b>TWS 32065</b>
black	<b>TWS 32066</b>

Labelling of identification labels



Item number
<b>TWS 32067</b>

Towel



Colour	Dimensions	Item number
blue	130 x 60 cm	<b>TWS 32068</b>
blue	140 x 100 cm	<b>TWS 32069</b>
blue	185 x 45 cm	<b>TWS 32062</b>
blue	185 x 60 cm	<b>TWS 32061</b>



**TWS 32070**

Quiver for instruments / HF-electrodes, 23 x 6 cm, polypropylene, autoclavable

Fixation clamp




Dimensions	Item number
klein/small	<b>TWS 32072</b>
groß/large	<b>TWS 32071</b>

Holding clamp




Ø	Item number
7 - 12 mm	<b>TWS 32076</b>
12 - 16 mm	<b>TWS 32077</b>
16 - 26 mm	<b>TWS 32073</b>
26 - 36 mm	<b>TWS 32074</b>
36 - 45 mm	<b>TWS 32075</b>

Holding pin for instruments to be used with wire baskets, without slot




Item number
<b>TWS 32078</b>

Holding pin for dividers to be used with wire baskets, with slot



Item number
<b>TWS 32079</b>

Dividers



Dimensions	Item number
50 x 20 mm	<b>TWS 32080</b>
130 x 20 mm	<b>TWS 32081</b>
225 x 20 mm	<b>TWS 32082</b>
460 x 20 mm	<b>TWS 32083</b>

Divider for sliding rails



Item number
<b>TWS 32084</b>

Sliding rail



Dimensions	Item number
230 x 26 mm	<b>TWS 32085</b>
470 x 26 mm	<b>TWS 32086</b>

Spiral holder



Dimensions

**126 x 30 mm**

**TWS 32087**

Spiral holder, double



Dimensions

**233 x 30 mm**

**TWS 32088**

Limiting bar



Dimensions

Item number

112 x 18 mm

**TWS 32089**

172 x 18 mm

**TWS 32090**

232 x 18 mm

**TWS 32091**

472 x 18 mm

**TWS 32092**

Limiting bar



Dimensions

Item number

230 x 26 mm

**TWS 32093**

470 x 26 mm

**TWS 32094**

Cleaning agent for anodized aluminium



Item number

**TWS 32095**

Foot for stainless steel trays, height 10 mm



Item number

**TWS 32096**

Foot for Mini and Special stainless steel trays, height 5 mm



Item number

**TWS 32097**

Silicone adhesive and caulk gun for spare gasket, 310 ml



Item number

**TWS 32098**

Spare silicone gasket for container lid



Container size	Length	Item number
1/1	158 cm	<b>TWS 32099</b>
3/4	136 cm	<b>TWS 32100</b>
1/2	101 cm	<b>TWS 32101</b>
Flach / Flat	101 cm	<b>TWS 32102</b>
Mini	77 cm	<b>TWS 32103</b>
Special	89 cm	<b>TWS 32104</b>

! NOT TO BE USED WITH DI-WATER! NOT REPROCESSABLE!



**TWS 32105 SET**

complete, with 8 tips and 3/4" coupling, 1.5 m pressure tube

**TWS 32106**

Water jet pistol, only



**TWS 32107**

for bottles



**TWS 32108**

for catheters



**TWS 32109**

for drainage tubes



**TWS 32110**

for receptacles



**TWS 32111**

for needles with Record-cone



**TWS 32112**

for needles with Luer



**TWS 32113**

for pipettes



**TWS 32114**

Cleaning tip "Water jet"



**TWS 32115**

Coupling 1/2"



**TWS 32116**

Coupling 3/4"

! SUITABLE FOR DI-WATER! REPROCESSABLE!



**TWS 32117 SET**

complete, with 4 tips, 1.5 m pressure tube

**TWS 32118**

Water jet pistol, only

Tips, SUITABLE FOR DI-WATER! REPROCESSABLE!



**TWS 32119**

Spray nozzle, stainless steel



**TWS 32120**

Luer-Lock connector



**TWS 32121**

Shower nozzle, stainless steel



**TWS 32122**

Luer-Lock connector, stainless steel



**TWS 32123**

Tube nozzle, size 6 (4 mm)



**TWS 32124**

Pressure tube for DI-Water, 1.5 m

Cleaning brushes with plastic handle and nylon brush



	Working length	Working length
Ø	330,0 mm	450,0 mm
3,0 mm	<b>TWS 32125</b>	<b>TWS 32128</b>
5,0 mm	<b>TWS 32126</b>	<b>TWS 32129</b>
10,0 mm	<b>TWS 32127</b>	<b>TWS 32130</b>

Cleaning brushes



	Total length
Ø	480,0 mm
3,0 mm	<b>TWS 32131</b>
5,0 mm	<b>TWS 32132</b>
7,0 mm	<b>TWS 32133</b>
12,0 mm	<b>TWS 32134</b>



**TWS 32135** 210 x 95 mm

Cleaning brush



**TWS 32136**

Nailbrush, autoclavable



**TWS 32137** 280 x 285 x 85 mm

Brush dispenser, for 20 brushes



**TWS 32138**

Protection tube for endoscopes up to Ø 5.0 mm and WL 188 mm



**TWS 32139**

Protection tube for endoscopes up to Ø 5.0 mm and WL 302 mm



KARL STORZ  
OLYMPUS

**TWS 32140**

Spare adapter for endoscopes



RICHARD WOLF

**TWS 32141**

Spare adapter for endoscopes



**TWS 32142**

Polish for endoscopes



**TWS 32143**

Protection cap for light entry of endoscopes



**TWS 32144**

Spare sealing cap

ULTRA-STOP Antifog solution

**TWS 32145**

30 ml, sterile, in pierce bottle

**TWS 32146**

25 ml, non sterile, in pipette bottle

**TWS 32147**

15 ml, non sterile, in spray bottle



**TWS 32148**

Luer-Lock male with hose coupling Ø 9 mm



**TWS 32149**

Luer-Lock cap, small




**TWS 32150**

Luer-Lock male with hose coupling Ø 6 mm




**TWS 32151**

Luer-Lock female with hose coupling Ø 6 mm



**TWS 32152**

Luer-Lock male with stopcock and hose coupling Ø 9 mm



**TWS 32153**

Luer-Lock male, rotating, with stopcock and Luer-Lock female



**TWS 32154**

Lubricant for stopcocks, 3 g