



Arm and Hand Surgery Table Finger Traps

TW Surgical Finger Traps: TW Surgical offers all sizes of Finger Traps to meet your surgical needs. To view additional Finger Trap sizes as follows:



Extra Small Finger Trap



Small Finger Trap



Medium Finger Trap



Large Finger Trap



Extra Large Finger Trap

TWS 1539/01	Extra Small Finger Trap Material: PE and Nylon
TWS 1539/02	Small Finger Trap Material: PE and Nylon
TWS 1539/03	Medium Finger Trap Material: PE and Nylon
TWS 1539/04	Large Finger Trap Material: PE and Nylon
TWS 1539/05	Extra Large Finger Trap Material: PE and Nylon



Arm and Hand Surgery Table

TW Surgical Finger Traps: TW Surgical offers all sizes of Finger Traps to meet your surgical needs. To view additional Finger Trap sizes as follows:




Product Description

-Stainless Steel Wire Mesh Finger Traction Device Set

(Sizes: XS, S, M, L)

-Re-usable Wire Finger Traps



TWS 1540/01	Extra Small Finger Trap Material: Stainless Steel
TWS 1540/02	Small Finger Trap Material: Stainless Steel
TWS 1540/03	Medium Finger Trap Material: Stainless Steel
TWS 1540/04	Large Finger Trap Material: Stainless Steel

Arm and Hand Surgery Table

Finger Traps



Triangle Traction Plate Device

Device can be used with IV poles or other extremity systems.

TWS 1525-01	Triangle Traction Plate Device Material: Aluminum and stainless steel
-------------	--------------------------------------------------------------------------

Horizontal Traction Tower, 18 Inches



Horizontal Traction Tower, 18 Inches, Solid Stainless Steel Construction

Adapts Easily to End of Arm Table, 3/4" Outside Diameter Square

Tube with 1/2" Telescoping Tube,

The strongest and most dependable, Adjustable height with clamp

Includes: Metal wire with metal loops

TWS 2335/01	Triangle Traction Tower 18 inches Material: Stainless steel
-------------	----------------------------------------------------------------

Arm and Hand Surgery Table

Finger Traction Assembly with Mounting Clamp

This finger traction assembly consists of durable stainless steel construction.

Use with 1-5 finger trap placements

Fully adjustable for both height and angulation requirements

Each placement wheel is independently secured in the desired position and adjust tension knob.

Available in a variety of sizes. Size.



TWS 2339/01

Triangle Traction Assembly

Material: Stainless steel