



## Micro Instrument – Micro Grip Forceps




TWS 11-1325/1205

Micro Grip Forceps 12cm Length, Straight, Fine Tips,  
0.5mm Thickness, 0.9mm Tip, Strong and Sharp Tip  
Stainless Steel and Silicone Rubber



TWS 11-1326/1257

Micro Grip Forceps 12.5cm Length, Straight, Round Fine Tips,  
0.07mm Thickness, 1.7mm Tip, Strong and Sharp Tip  
Stainless Steel and Silicone Rubber



TWS 11-1327/1211


Micro Grip Forceps 12cm Length, Straight, Round Fine Tips,  
0.11mm Thickness, 1.5mm Tip, Strong and Sharp  
Stainless Steel and Silicone Rubber

**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom.

info@twrightmedical.com


## Micro Instrument – Micro Grip Forceps



TWS 11-1328/1259	Micro Grip Forceps 12.5cm Length, Straight, Extra Fine Tips, 0.09mm Thickness, 0.16mm Tip, Strong and Sharp Stainless Steel and Silicone Rubber
TWS 11-1328/1109	Micro Grip Forceps 11cm Length, Straight, Extra Fine Tips, 0.09mm Thickness, 0.2mm Tip, Strong and Sharp Stainless Steel and Silicone Rubber



TWS 11-1329/1303	Micro Grip Forceps 13cm Length, Straight, Fine Tips, 0.3mm Thickness, 0.25mm Tips, Strong and Sharp Stainless Steel and Silicone Rubber
------------------	--




TWS 11-1330/1109	Micro Grip Forceps 11cm Length, Straight, Extra Fine Tips, 0.09mm Thickness, 0.11mm Tip, Strong and Sharp Stainless Steel and Silicone Rubber
------------------	--


**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom.  
info@twrightmedical.com


## Micro Instrument – Micro Grip Forceps



TWS 11-1331/1106	Micro Grip Forceps 11cm Length, Straight, Extra Fine Tips, 0.06mm Thickness, 0.1mm Tip, Strong and Sharp Stainless Steel and Silicone Rubber
------------------	--




TWS 11-1332/1158	Micro Grip Forceps 11.5cm Length, Curved, Extra Fine Tips, 0.08mm Thickness, 0.12mm Bent Tip, Strong and Sharp Stainless Steel and Silicone Rubber
------------------	--




TWS 11-1333/1151	Micro Grip Forceps 11.5cm Length, Curved, Extra Fine Tips, 0.1mm Thickness, 0.15mm Bent Tip, Strong and Sharp Stainless Steel and Silicone Rubber
------------------	---

## Micro Instrument – Micro Grip Forceps




TWS 11-1334/1154

Micro Grip Forceps 11.5cm Length, Curved, Extra Fine Tips, 0.04mm Thickness, 0.12mm Bent Tip, Strong and Sharp Stainless Steel and Silicone Rubber



TWS 11-1335/1209

Micro Grip Forceps 12cm Length, Curved, Extra Fine Tips, 0.09mm Thickness, 0.16mm Angled Tip, Strong and Sharp Stainless Steel and Silicone Rubber



TWS 11-1336/1155


Micro Grip Forceps 11.5cm Length, Curved, Fine Tips, 4.6mm Thickness, 11mm Angled Blade 105°, Strong and Sharp Stainless Steel and Silicone Rubber

**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom.


info@twrightmedical.com

## Micro Instrument – Micro Grip Forceps




TWS 11-1337/1155

Micro Grip Forceps 11.5cm Length, Curved, Fine Tips, 4.2mm Thickness, 6mm Angled Blades 105°, Strong and Sharp Stainless Steel and Silicone Rubber



TWS 11-1338/1285

Micro Grip Forceps 12cm Length, Curved, Fine Tips, 8.5mm Thickness, 12mm Angled Blades 135°, Strong and Sharp Stainless Steel and Silicone Rubber



TWS 11-1339/1151

Micro Grip Forceps 11.5cm Length, Curved, Fine Tips, 9.5mm Thickness, 10mm Angled Blades 135°, Strong and Sharp Stainless Steel and Silicone Rubber

**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom.

info@twrightmedical.com