


Collis Posterior Lumbar Retractor Set



Collis Retractor Frame and Universal Swivel Case 1


The set includes 32 pairs of blades

32 available sizes

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set



Collis Retractor Frame and Universal Swivel Case 2

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set




Collis Retractor Frame and Universal Swivel Case 3

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set

The Collis Retractor Frame and Universal Swivel is a component of a surgical system used primarily in **posterior lumbar spinal surgery** to hold back muscles and ligaments, providing surgeons with clear, unobstructed access to the spinal column for procedures like disc removal or fusion, with the frame lying flat to minimize torque and the swivel allowing for precise blade positioning.



Functions:

- **Tissue Retraction:** It holds the paraspinal muscles and ligaments away from the surgical site.
- **Enhanced Visualization:** Creates a clear operative field for surgeons to see and operate on the spine.
- **Stability:** The frame's low-profile design helps prevent torque on the spine.
- **Versatility:** Used in various spinal surgeries, from microdiscectomies to scoliosis corrections.

Surgical Applications

The Collis Retractor system is used for a variety of posterior lumbar spinal procedures, accommodating different incision sizes and depths.

- **Microsurgery:** It is suitable for delicate procedures such as those for unilateral disc protrusion.
- **Decompression and Fusion:** It is used in larger operations, including decompression and fusion for conditions like scoliosis.
- **Hemilaminectomy:** The set includes specific blades designed for use in hemilaminectomy procedures.

Features of the System

- **Self-retaining:** The system is designed to hold tissues back without constant manual force, freeing up the surgeon's hands and providing stable, consistent exposure.
- **Versatile Components:** The frame works with a large selection of interchangeable blades (available in eight lengths and four widths) to adapt to different patient anatomies.
- **Enhanced Visualization:** By holding the retractor flat on the patient's back, it ensures an unobstructed view of the operative field, which can be further improved with optional light guides.
- **Reduced Trauma:** The design aims to minimize tissue trauma by providing stable, adjustable retraction that avoids excessive pressure on any single area when used correctly.

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com


Collis Posterior Lumbar Retractor Set

Collis Retractor Frame and Universal Swivel serve as the foundational, stable component to which various blades are attached to perform these specific spinal surgeries.

The Collis Posterior Lumbar Retractor Blades Set offers a large selection of blades in eight lengths and four widths, for retracting the paraspinal muscles and ligaments on both sides of the spinal column during any type of spinal surgery.

The set includes 32 pairs of blades that cover the gamut of procedure sizes: from microsurgery for unilateral disc protrusion, to large decompression for fusion of scoliosis.

Working with 32 available sizes ensures having the right blade for any patient or opening size, so the retractor lies flat on the patient's back for an unobstructed view of the operative field. Torque is reduced as well.



TWS 25-3000-01 Sterilizing Case

TWS 25-2020-1000 Collis Posterior Lumbar Retractor frame

Collis Retractor Frame and Universal Swivel Case 1

TWS 25-2020-0000	Collis Retractor Frame and Universal Swivel System
TWS 25-2020-1000	Collis Posterior Lumbar Retractor frame with 6 nuts Material: Stainless Steel
TWS 25-2020-100	Universal swivel for 5.5 mm post Material: Stainless Steel
TWS 25-3000-01	Sterilizing Case for Posterior Lumbar Retractor, Stainless Steel

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set



Collis Retractor Hemilaminectomy Blade



Curved Nerve Root Retractor Blades


Collis Retractor Frame and Universal Swivel Case 1

TWS 25-2030-01	Hemilaminectomy Blade; 1" (2.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-02	Hemilaminectomy Blade; 1 3/8" (3.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-03	Hemilaminectomy Blade; 1 3/4" (4.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-04	Hemilaminectomy Blade; 2 1/8" (5.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-05	Hemilaminectomy Blade; 2 1/2" (6.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-06	Hemilaminectomy Blade; 3" (7.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-07	Hemilaminectomy Blade; 3 3/8" (8.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2030-08	Hemilaminectomy Blade; 3 3/4" (9.5 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2050-01	Curved Nerve Root Retractor; 7" (18 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2050-02	Curved Nerve Root Retractor; 9" (23 cm) Length x 1/2" (1.2 cm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2050-03	Curved Nerve Root Retractor; 7" (18 cm) Length x 5/16" (8 mm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel
TWS 25-2050-04	Curved Nerve Root Retractor; 9" (23 cm) Length x 5/16" (8 mm) Wide Collis Posterior Lumbar Retractor Blade, Stainless Steel

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set




TWS 25-3000-02 Sterilizing Case

Collis Retractor Frame and Universal Swivel Case 2



Collis Retractor Blades



Collis Retractor Blades

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set



Collis Retractor Blades Case 2

TWS 25-2020-01	Collis Retractor Blades Size 1, Posterior Lumbar Retractor, Stainless Steel 1 1/2" (4 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-02	Collis Retractor Blades Size 2, Posterior Lumbar Retractor, Stainless Steel 2" (5 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-03	Collis Retractor Blades Size 3, Posterior Lumbar Retractor, Stainless Steel 2 3/8" (6 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-04	Collis Retractor Blades Size 4, Posterior Lumbar Retractor, Stainless Steel 2 3/4" (7 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-05	Collis Retractor Blades Size 5, Posterior Lumbar Retractor, Stainless Steel 3 1/8" (8 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-06	Collis Retractor Blades Size 6, Posterior Lumbar Retractor, Stainless Steel 3 1/2" (9 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-07	Collis Retractor Blades Size 7, Posterior Lumbar Retractor, Stainless Steel 4" (10 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-08	Collis Retractor Blades Size 8, Posterior Lumbar Retractor, Stainless Steel 4 3/8" (11 cm) Length x 1/2" (1.2 cm) Wide x 2 Nos.
TWS 25-2020-09	Collis Retractor Blades Size 9, Posterior Lumbar Retractor, Stainless Steel 1 1/2" (4 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-10	Collis Retractor Blades Size 10. Posterior Lumbar Retractor, Stainless Steel 2" (5 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-11	Collis Retractor Blades Size 11, Posterior Lumbar Retractor, Stainless Steel 2 3/8" (6 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-12	Collis Retractor Blades Size 12, Posterior Lumbar Retractor, Stainless Steel 2 3/4" (7 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-13	Collis Retractor Blades Size 13, Posterior Lumbar Retractor, Stainless Steel 3 1/8" (8 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-14	Collis Retractor Blades Size 14, Posterior Lumbar Retractor, Stainless Steel 3 1/2" (9 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-15	Collis Retractor Blades Size 15, Posterior Lumbar Retractor, Stainless Steel 4" (10 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-2020-16	Collis Retractor Blades Size 16, Posterior Lumbar Retractor, Stainless Steel 4 3/8" (11 cm) Length x 1" (2.5 cm) Wide x 2 Nos.
TWS 25-3000-02	Sterilizing Case for Posterior Lumbar Retractor, Stainless Steel

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com

Collis Posterior Lumbar Retractor Set




TWS 25-3000-03 Sterilizing Case

Collis Retractor Blades Case 3



Collis Retractor Blades



Collis Retractor Blades

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.
info@twrightmedical.com



TW Surgical

Collis Posterior Lumbar Retractor Set



Collis Retractor Blades Case 3

TWS 25-2020-17	Collis Retractor Blades Size 17, Posterior Lumbar Retractor, Stainless Steel 1 1/2" (4 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-18	Collis Retractor Blades Size 18, Posterior Lumbar Retractor, Stainless Steel 2" (5 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-19	Collis Retractor Blades Size 19, Posterior Lumbar Retractor, Stainless Steel 2 3/8" (6 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-20	Collis Retractor Blades Size 20, Posterior Lumbar Retractor, Stainless Steel 2 3/4" (7 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-21	Collis Retractor Blades Size 21, Posterior Lumbar Retractor, Stainless Steel 3 1/8" (8 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-22	Collis Retractor Blades Size 22, Posterior Lumbar Retractor, Stainless Steel 3 1/2" (9 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-23	Collis Retractor Blades Size 23, Posterior Lumbar Retractor, Stainless Steel 4" (10 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-24	Collis Retractor Blades Size 24, Posterior Lumbar Retractor, Stainless Steel 4 3/8" (11 cm) Length x 1 1/2" (3.8 cm) Wide x 2 Nos.
TWS 25-2020-25	Collis Retractor Blades Size 25, Posterior Lumbar Retractor, Stainless Steel 1 1/2" (4 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-26	Collis Retractor Blades Size 26, Posterior Lumbar Retractor, Stainless Steel 2" (5 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-27	Collis Retractor Blades Size 27, Posterior Lumbar Retractor, Stainless Steel 2 3/8" (6 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-28	Collis Retractor Blades Size 28, Posterior Lumbar Retractor, Stainless Steel 2 3/4" (7 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-29	Collis Retractor Blades Size 29, Posterior Lumbar Retractor, Stainless Steel 3 1/8" (8 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-30	Collis Retractor Blades Size 30, Posterior Lumbar Retractor, Stainless Steel 3 1/2" (9 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-31	Collis Retractor Blades Size 31, Posterior Lumbar Retractor, Stainless Steel 4" (10 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-2020-32	Collis Retractor Blades Size 32, Posterior Lumbar Retractor, Stainless Steel 4 3/8" (11 cm) Length x 2" (5 cm) Wide x 2 Nos.
TWS 25-3000-03	Sterilizing Case for Posterior Lumbar Retractor, Stainless Steel

Treasure Wright (U.K.) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.

info@twrightmedical.com