

DESCRIPTION

The Taylor retractor is designed with a single sharp tip at the end of the blade to let the intervertebral disc rongeurs to nestle down in it. The Taylor Small Blades provide a low-profile surgical retraction and allow placing the retractor just on the skin without any space between the retractor and skin by sliding the retractor over its blade downwards.



TAYLOR RETRACTOR				
CAT. NO.	DESCRIPTION	BLADE WIDTH	SIZE	
			BLADE LENGTH	
TWS 74-2519/50	TAYLOR RETRACTOR, SMALL	19mm	50mm	
TWS 74-2519/75	TAYLOR RETRACTOR, SMALL	19mm	75mm	
TWS 74-2519/100	TAYLOR RETRACTOR, SMALL	19mm	100mm	
TWS 74-2519/125	TAYLOR RETRACTOR, SMALL	19mm	125mm	
TWS 74-2519/150	TAYLOR RETRACTOR, SMALL	19mm	150mm	



TAYLOR RETRACTOR				
CAT. NO.	DESCRIPTION	BLADE WIDTH	SIZE	
			BLADE LENGTH	
TWS 74-2520/75	TAYLOR RETRACTOR, STANDARD	30mm	75mm	
TWS 74-2520/100	TAYLOR RETRACTOR, STANDARD	30mm	100mm	
TWS 74-2520/150	TAYLOR RETRACTOR, STANDARD	30mm	150mm	



TAYLOR RETRACTOR				
CAT. NO.	DESCRIPTION	BLADE WIDTH	SIZE	
			BLADE LENGTH	
TWS 74-2530/75	TAYLOR RETRACTOR, SMALL, CONCAVE	19mm	75mm	
TWS 74-2530/100	TAYLOR RETRACTOR, SMALL, CONCAVE	19mm	100mm	



TAYLOR RETRACTOR				
CAT. NO.	DESCRIPTION	BLADE WIDTH	SIZE	
			BLADE LENGTH	
TWS 74-2540/75	TAYLOR RETRACTOR, SMALL, FIBEROPTIC	19mm	75mm	
TWS 74-2540/100	TAYLOR RETRACTOR, SMALL, FIBEROPTIC	19mm	100mm	