

## Jackson-Burrows self-retaining retractor

The Jackson-Burrows self-retaining retractor is a specialized surgical instrument, commonly used in orthopaedic and general surgeries to hold back tissues, skin, or organs for enhanced visibility. A self-locking mechanism and available with different blade depths providing access to deep surgical fields.

Jackson Burrows Retractors Blunt is a finger ring, self-retaining retractor used during intracranial procedures to hold back the cerebellum. It is available with blunt interlocking prongs in various configurations, angled or straight arms, and multiple lengths.



- **Purpose:** Primarily designed for retracting wound and incision edges to provide optimal surgical view.
- **Design:** Self-retaining, allowing for hands-free, secure tissue retraction.
- **Structure:** Sharp and Blunt hooked jaws or prongs and an ergonomic handle
- **Material:** Typically crafted from high-grade stainless steel for durability.
- **Blunt/Sharp:** Available with blunt jaws for delicate tissues or sharp prongs.
- **Size:** Available with different deep sizes for deep or narrow surgical fields.

**Mechanism:** Features a ratchet lock that allows surgeons to secure the instrument in place, maintaining tension on the wound edges hands-free.

**Design:** Typically includes 1 x 1 sharp and blunt prongs (one prong on each arm) designed to minimize trauma to delicate tissues while providing secure retraction.

**Treasure Wright (U.K.) Ltd.**

128-160 City Road, London, EC1V 2NX, United Kingdom.

info@twrightmedical.com

## Jackson-Burrows self-retaining retractor



TWS 160-1750/10S	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 20mm Jaw Depth, 175mm Long, Sharp Jaw Tip 10mm length Material: Stainless Steel
TWS 160-1751/10S	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 30mm Jaw Depth, 175mm Long, Sharp Jaw Tip 10mm length Material: Stainless Steel
TWS 160-1752/10S	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 40mm Jaw Depth, 175mm Long, Sharp Jaw Tip 10mm length Material: Stainless Steel
TWS 160-1750/10B	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 20mm Jaw Depth, 175mm Long, Blunt Jaw Tip 10mm length Material: Stainless Steel
TWS 160-1751/10B	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 30mm Jaw Depth, 175mm Long, Blunt Jaw Tip 10mm length Material: Stainless Steel
TWS 160-1752/10B	Jackson-Burrows self-retaining retractor with Ratchet Lock 1 x 1 prong, 40mm Jaw Depth, 175mm Long, Blunt Jaw Tip 10mm length Material: Stainless Steel

### Treasure Wright (U.K.) Ltd.