

Gel Positioners




TWS 602075/75



TWS 612075/75



TWS 622075/45



TWS 631919/06

CLOSED HEAD RING

Code	LxWxH cm	LxWxH inch
TWS 60904/02	9 x 4 x 2	3.5 x 1.6 x 0.8
TWS 601406/35	14 x 6 x 3.5	5.5 x 2.4 x 1.4
TWS 602075/45	20 x 7.5 x 4.5	7.9 x 3.0 x 1.8
TWS 602075/75	20 x 7.5 x 7.5	7.9 x 3.0 x 3.0

OPEN HEAD RING

Code	LxWxH cm	LxWxH inch
TWS 610904/02	9 x 4 x 2	3.5 x 1.6 x 0.8
TWS 611406/35	14 x 6 x 3.5	5.5 x 2.4 x 1.4
TWS 612075/45	20 x 7.5 x 4.5	7.9 x 3.0 x 1.8
TWS 612075/75	20 x 7.5 x 7.5	7.9 x 3.0 x 3.0

PADLESS OPEN HEAD RING

Code	LxWxH cm	LxWxH inch
TWS 620904/02	9 x 4 x 2	3.5 x 1.6 x 0.8
TWS 621406/35	14 x 6 x 3.5	5.5 x 2.4 x 1.4
TWS 622075/45	20 x 7.5 x 4.5	7.9 x 3.0 x 1.8

CONTOURED HEADREST


Code	LxWxH cm	LxWxH inch
TWS 631919/06	19 x 19 x 6	7.5 x 7.5 x 2.4




TWS 642319/13



TWS 652723/07



TWS 662824/09



TWS 671606/16

Code	LxWxH cm	LxWxH inch
TWS 642823/14	28 x 23 x 14	11.0 x 9.1 x 5.7
TWS 642319/13	23 x 19 x 13	9.1 x 7.5 x 5.1

SUPINE HEAD REST

Code	LxWxH cm	LxWxH inch
TWS 662824/09	28 x 24 x 9	11 x 9.4 x 3.5

OPHTHALMIC HEADREST

Code	LxWxH cm	LxWxH inch
TWS 652723/07	27 x 23 x 7	10.6 x 9.1 x 2.8

FACE PAD


Code	LxWxH cm	LxWxH inch
TWS 671606/16	16 x 6 x 1.6	6.3 x 2.4 x 0.6

Treasure Wright (UK) Ltd.

128-160 City Road, London, EC1V 2NX, United Kingdom.

Email: info@twrightmedical.com

Gel Positioners



ARMBOARD PAD

Useful during various surgical procedures in all positions, these pads provide both comfort and excellent pressure relief.

CONTOURED ARMBOARD PAD WITH FOAM BASE

The contour design on these arm board pads provide superior stability and support.

ARM PROTECTOR

Mainly used in the lateral position, the arm protector comes with adjustable Velcro to provide stable pressure relief.

FASTENING BELT

The Gel Fastening Straps hold the patients arms on the table pad. The easily adjustable straps are made of soft polyurethane gel that protects elbows and ankles from pressure.



TWS 106012/02

TWS 116315/06

TWS 124621/08

TWS 133607/15

ARMBOARD PAD

ARM PROTECTOR

CONTOURED ARMBOARD PAD

FASTENING BELT

Code	LxWxH cm	LxWxH inch
TWS 104012/02	40 x 12.5 x 2	15.7 x 4.9 x 0.8
TWS 106012/02	60 x 12.5 x 2	23.6 x 4.9 x 0.8

Code	LxWxH cm	LxWxH inch
TWS 114015/06	40 x 15 x 6	15.7 x 5.9 x 2.4
TWS 115015/06	50 x 15 x 6	19.7 x 5.9 x 2.4
TWS 116015/06	60 x 15 x 6	23.6 x 5.9 x 2.4
TWS 116315/06	63 x 15 x 6	24.8 x 5.9 x 2.4

Code	LxWxH cm	LxWxH inch
TWS 123812/08	38 x 12 x 0.8	15.0 x 4.7 x 0.3
TWS 124621/08	46 x 21 x 0.8	18.1 x 8.3 x 0.3

Code	LxWxH cm	LxWxH inch
TWS 133307/15	33 x 7 x 1.5	13.0 x 2.8 x 0.6
TWS 133607/15	36 x 7 x 1.5	14.2 x 2.8 x 0.6

Gel Positioners



TWS 145018/15



TWS 153007/75



TWS 165016/13



TWS 175012/07

CHEST ROLL (ADULT)

Code	LxWxH cm	LxWxH inch
TWS 143010/85	30 x 10 x 8.5	11.8 x 3.9 x 3.3
TWS 141010/85	40 x 10 x 8.5	15.7 x 3.9 x 3.3
TWS 145010/85	50 x 10 x 8.5	19.7 x 3.9 x 3.3
TWS 143015/10	30 x 15 x 10	11.8 x 5.9 x 3.9
TWS 144015/10	40 x 15 x 10	15.7 x 5.9 x 3.9
TWS 145015/10	50 x 15 x 10	19.7 x 5.9 x 3.9
TWS 145018/15	50 x 18 x 15	19.7 x 7.1 x 5.9

CHEST ROLL (PEDIATRIC)

Code	LxWxH cm	LxWxH inch
TWS 151505/05	15 x 5 x 5	5.9 x 2.0 x 2.0
TWS 152005/05	20 x 5 x 5	7.9 x 2.0 x 2.0
TWS 152505/05	25 x 5 x 5	9.8 x 2.0 x 2.0
TWS 153005/05	30 x 5 x 5	11.8 x 2.0 x 2.0
TWS 151507/75	15 x 7.5 x 7.5	5.9 x 3.0 x 3.0
TWS 152007/75	20 x 7.5 x 7.5	7.9 x 3.0 x 3.0
TWS 152507/75	25 x 7.5 x 7.5	9.8 x 3.0 x 3.0
TWS 153007/75	30 x 7.5 x 7.5	11.8 x 3.0 x 3.0

CONTOURED CHEST ROLL

Code	LxWxH cm	LxWxH inch
TWS 163016/13	30 x 16 x 13	11.8 x 6.3 x 5.1
TWS 164016/13	40 x 16 x 13	15.7 x 6.3 x 5.1
TWS 165016/13	50 x 16 x 13	19.7 x 6.3 x 5.1

UNIVERSAL PATIENT POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 171608/45	16 x 8 x 4.5	6.3 x 3.1 x 1.8
TWS 173305/05	33 x 7 x 5	13.0 x 2.8 x 2.0
TWS 173312/06	33 x 12 x 6	13.0 x 4.7 x 2.4
TWS 173612/07	36 x 12 x 7	14.2 x 4.9 x 2.8
TWS 172512/05	25 x 12 x 5	9.8 x 4.9 x 2.0
TWS 173812/06	38 x 12 x 6	15.0 x 4.9 x 2.5
TWS 175012/07	50 x 12 x 7	19.7 x 4.9 x 2.8

Gel Positioners



TWS 185515/06R L



TWS 197150/15



TWS 206015/10



TWS 215085/08

CHEST POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 187120/75R	71 x 20 x 7.5	28.0 x 7.9 x 3.0
TWS 187120/75L	71 x 20 x 7.6	28.0 x 7.9 x 3.0
TWS 185515/06R	55 x 15 x 6	21.6 x 3.1 x 1.2
TWS 185515/06L	55 x 15 x 6	21.6 x 3.1 x 1.2

PATIENT POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 201608/06	16 x 8 x 6	6.3 x 3.1 x 2.4 (2/gel, 4/foam)
TWS 207020/08	70 x 20 x 8	27.6 x 7.9 x 3.1 (2/gel, 6/foam)
TWS 203015/10	30 x 15 x 10	11.8 x 5.9 x 3.9 (2/gel, 8/foam)
TWS 205015/10	50 x 15 x 10	19.7 x 5.9 x 3.9 (2/gel, 8/foam)
TWS 206015/10	60 x 15 x 10	23.6 x 5.9 x 3.9 (2/gel, 8/foam)

LATERAL POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 197150/15	71 x 50 x 15.2	28.0 x 19.7 x 6.0

PATIENT POSITIONING BELT

Code	LxWxH cm	LxWxH inch
TWS 215085/08	50 x 8.5 x 0.8	19.9 x 3.3 x 0.3

Gel Positioners



TWS 222011/07



TWS 231710/05



TWS 241718/13



TWS 251723/13

HEEL SUPPORT

Code	LxWxH cm	LxWxH inch
TWS 222011/07	20 x 11 x 7	7.9 x 4.3 x 2.8

HEEL CUP

Code	LxWxH cm	LxWxH inch
TWS 231710/05	17 x 10 x 5	6.7 x 4.1 x 2.0

COMFORT HEEL CUP

Code	LxWxH cm	LxWxH inch
TWS 241718/13	17 x 18 x 1.3	6.7 x 7.1 x 0.5

COMFORT HEEL PAD

Code	LxWxH cm	LxWxH inch
TWS 251723/13	17 x 23 x 1.3	6.7 x 9.1 x 0.5

Gel Positioners




TWS 405552/01



TWS 412328/13



TWS 423623/38



TWS 433724/05

CRUTCH STIRRUP PAD

Code	LxWxH cm	LxWxH inch
TWS 403123/01	31 x 23 x 1	12.4 x 9.3 x 0.4
TWS 403924/01	39 x 24 x 1	15.6 x 9.6 x 0.4
TWS 405230/01	52 x 30 x 1	20.5 x 11.8 x 0.4
TWS 405552/01	55 x 52 x 1	21.7 x 20.5 x 0.4

BOOT STIRRUP PAD

Code	LxWxH cm	LxWxH inch
TWS 412328/13	23 x 28 x 1.3	9.1 x 11.0 x 0.5

GRADUATED THORACIC POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 423623/38	36 x 23 x 3.8	14.2 x 9.1 x 1.5

TRACTION BOOT PAD

Code	LxWxH cm	LxWxH inch
TWS 433724/05	37 x 24 x 0.5	14.6 x 9.4 x 0.2

Gel Positioners



TWS 443515/75

DOME POSITIONER

Code	LxWxH cm	LxWxH inch
TWS 443515/75	35 x 15 x 7.5	13.8 x 5.9 x 3



TWS 454445/01



TWS 467112/08

ANKLE STRAP


Code	LxWxH cm	LxWxH inch
TWS 453104/01	31 x 4 x 1	12.2 x 1.6 x 0.4
TWS 454445/01	44 x 4.5 x 1	17.3 x 1.8 x 0.4

LITHOTOMY POLE COVER

Code	LxWxH cm	LxWxH inch
TWS 467112/08	71 x 12 x 0.8	28.0 x 4.7 x 0.3

Gel Positioners

CUSHION



Gel cushion

The gel cushion has a smooth surface, effectively dissipating heat to provide excellent pressure distribution. The high-quality polyurethane gel is durable and allows for easy cleaning and upkeep.

Gel cushion

Code	LxWxH cm	LxWxH inch
TWS 504141/16	41x41x1.6	16x16x0.63
TWS 504641/16	46x41x1.6	18x16x0.63
TWS 504646/16	46x46x1.6	18x18x0.63
TWS 505146/16	51x46x1.6	20x18x0.63

Gel & Form cushion 1

Foam Cushion v1 provides softness and comfortability, with its plate-shaped foam covered with a layer of pressure-relieving polyurethane gel. The foam is 4cm thick, giving adequate support to patients at risk.

Gel & Form cushion 1

Code	LxWxH cm	LxWxH inch
TWS 504141/56	41x41x5.6	16x16x2.2
TWS 504641/56	46x41x5.6	18x16x2.2
TWS 504646/56	46x46x5.6	18x18x2.2
TWS 505146/56	51x46x5.6	20x18x2.2

Gel Positioners



Gel & Form cushion 2

Foam Cushion 2 is a two-dimensional foam with a hollowed-out design to provide better comfort to user's thighs. Pressure in the wider side is equally distributed, making this cushion ideal for user's who spend a long time sitting in the same place.

Gel & Form cushion 2

Code	L×W×H cm	L×W×H inch
TWS 514141/07	41×41×7	16×16×2.8
TWS 514641/07	46×41×7	18×16×2.8
TWS 514646/07	46×46×7	18×18×2.8
TWS 515146/07	51×46×7	20×18×2.8

Gel & Form cushion 3

Foam Cushion 3 is a three-dimensional foam, designed not only to hug the thighs but the hips as well. This outstanding feature provides maximum pressure distribution and bedsores protection in the hips. High risk patients are advised to use this to cushion their wheelchairs.

Gel & Form cushion 3

Code	L×W×H cm	L×W×H inch
TWS 514141/09	41×41×9	16×16×3.5
TWS 514641/09	46×41×9	18×16×3.5
TWS 514646/09	46×46×9	18×18×3.5
TWS 515146/09	51×46×9	20×18×3.5